

C3 at 

Complement C3

Turbidimetric method

<p>REF 3170005 1 x 50 mL</p> <p>CONTENTS R1.Reagent 1 x 50 mL</p>
<p>For <i>in vitro</i> diagnostic use only</p>

PRINCIPLE

C3 at is a quantitative turbidimetric assay for the measurement of the component complement C3 in human serum or plasma.

Anti-human C3 antibodies form insoluble complexes when mixed with samples containing C3. The scattering light of the immunocomplexes depends of the C3 concentration in the patient sample, and can be quantified by comparison from a calibrator of known C3 concentration.

REAGENTS COMPOSITION

R1 **C3 at.** Goat antibodies anti-human C3, tris buffer 20 mmol/L, pH 8.2. Sodium azide 0.95 g/L.

Plasma Protein Multicalibrator. Protein Calibrator. Optional . Ref: 3910005.

Precautions: The reagents contain sodium azide 0.95 g/L. Avoid any contact with skin or mucous.

STORAGE AND STABILITY

- Store at 2-8°C.
The reagent is stable until the expiry date stated on the label, when stored tightly closed at 2-8°C and contaminations are prevented during their use. Does not use the reagent after the expiry date.
- Presence of particles, turbidity and/or the absorbance of blank reagent > 0.3 at 340 nm are sign of deterioration.

REAGENT PREPARATION

R1 Ready to use.

Calibration curve. Dilute the Plasma Protein Calibrator in NaCl 9 g/L as follow:

Dilution	1	2	3	4	5	6
Calibrator (µL)	--	10	25	50	75	100
NaCl 9 g/L (µL)	100	90	75	50	25	--
Factor	0	0.1	0.25	0.5	0.75	1.0

Multiply the concentration of the C3 Protein Calibrator by the corresponding factor to obtain the C3 concentration of each dilution.

SAMPLES

Fresh serum and EDTA or heparinized plasma. C3 in serum or plasma is stable 7 days at 2-8°C or 3 months at -20°C.

Samples with presence of fibrin should be centrifuged before testing. Highly hemolyzed or lipemic samples are not suitable for testing.

INTERFERENCES

Bilirubin (10 mg/dL), hemoglobin (4 g/L) and rheumatoid factors (300 UI/mL) may affect the results. Lipemia (12 g/L) does not interfere. Other substances may interfere⁵.

MATERIAL REQUIRED

- Thermostatic bath at 37°C.
- Spectrophotometer or photometer thermostatable at 37°C capable to read at 340 ± 20 nm.
- Cuvettes with 1cm pathlength.
- Pipettes to measure reagent and samples.

PROCEDURE

- Prewarm the reagent and the photometer (cuvette holder) to 37°C.
- Using distilled water zero the instrument at 340 nm.
- Pipette into a cuvette:

Sample / Calibrator	10 µL
Reagent (R1)	1.0 mL

- Mix well and insert the cuvette into the photometer. Record the absorbance (A) after 2 minutes of the sample or calibrator addition.

CALCULATION

Plot the different absorbance values (A) against the C3 concentration of each calibrator dilution. C3 concentration in the sample is calculated by interpolation of its (A) value in the calibration curve.



REFERENCE VALUES

Adults³: 90 – 180 mg/dL

Newborn⁴: 70 – 196 mg/dL

It is recommended that each laboratory establishes its own reference range.

QUALITY CONTROLS

To ensure adequate quality control (QC), each run should include a set of controls (normal and abnormal) with assayed values handled as unknowns.

REF 3915010 PLASMA PROTEIN CONTROL N-I
Normal level. Assayed.

REF 3915015 PLASMA PROTEIN CONTROL N-II
Abnormal level. Assayed.

Each laboratory should establish its own Quality Control scheme and corrective actions if controls do not meet the acceptable tolerances.

CLINICAL SIGNIFICANCE

C3 is the complement component of highest concentration in plasma and the most functionally linked between the classical and alternative pathway activations¹. Hepatic cells synthesize C3, although bacterial endotoxins induce synthesis through monocytes and fibroblasts.

Increased and decreased levels of C3 both have clinical significance.

Increased levels are closely related with acute-phase response (trauma, inflammatory process), biliary obstruction and focal glomeruloesclerosis⁴.

Decreased levels are related with genetic deficiency (risk for infection, particularly by encapsulated bacteria), or acquired deficiency (collagen vascular diseases and severe infections)⁴.

ANALYTICAL PERFORMANCE

- **Linearity limit.** Up to 500 mg/dL, under the described assay conditions. Samples with higher concentrations should be diluted 1/5 en CINA 9 g/L and retested again.
- **Detection limit.** Values less than 4.2 mg/dL give non-reproducible results.
- **Analytical sensitivity.** Using this reagent and method an ΔA of 3.42 mA at 340 nm is equivalent to 1 mg/dL of C3 at a concentration of 192 mg/dL.
- **Prozone effect.** Prozone effect is not observed up to 800 mg/dL.

- Precision.

mg/dL	Within-run		Between-run	
Mean	79.3	151.2	79.3	151.2
SD	3.2	6.5	3.8	10.5
CV%	4.0	4.2	4.8	6.9
N	10	10	10	10

Instrument: Cobas Mira

- **Accuracy:** Results obtained with this reagent did not show systematic differences when compared with commercial reagents of similar characteristics. Details of comparison are available on request.

NOTES

1. This method may be used with different instruments. Any application to an instrument should be validated to demonstrate that results meet the performance characteristics of the method. It is recommended to validate periodically the instrument. Contact to the distributor for any question on the application method.
2. The linearity limit depends on the sample/reagent ratio, as well as the analyzer used. It will be higher by decreasing the sample volume, although the sensitivity of the test will be proportionally decreased.
3. Clinical diagnosis should not be made on findings of a single test result, but should integrate both clinical and laboratory data.

REFERENCES

1. Lambris John D. *Immunology Today* 9: 387-393 (1988).
2. Price CP et al. *Ann Clin Biochem* 20: 1-14 (1983).
3. Dati F et al. *Eur J Clin Chem Clin Biochem* 34 : 517-520 (1996).
4. Tietz *Textbook of Clinical Chemistry*, 3rd Ed. Burtis CA, Ashwood ER. WB Saunders Co., (1999).
5. Young DS. *Effects of drugs on clinical laboratory tests*. 3th ed. AACC Press (1997).
6. Friedman and Young. *Effects of the disease on clinical laboratory tests*, 3th ed. AACC Press, 1997.

