IgA + IgG + IgM id 

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IgA + IgG + IgM

Identification of the protein by radial immunodiffusion

PRINCIPLE

The examined protein, diffusing in agarose gel containing a specific antibody will form an immuno-complex, visible as a ring around the well. The ring diameter is in proportion to the concentration of the analyzed protein. The proportion corresponds to the diffusion time. In fact, at the end (72 h for IgA and IgG, 96 h for IgM), the square of the diameter will be in linear proportion to the concentration (*procedure 1*), while after a shorter period of time the square of the diameter will be in a logarithmic proportion to the concentration (*procedure 2*). In both cases, a calibration curve should be constructed, using at least three calibration points. However, a reference table is provided showing the relation between any concentration and the diameter at the end of the analysis time.

REAGENT COMPOSITION

Plate with agarose gel containing the antiserum for the specific proteins to be analyzed.

STORAGE AND STABILITY

Store at 2-8°C.
The plates are ready to be used.
The Reagents are stable until the expiry date stated on the label.

SAMPLES

Serum or plasma. Stable six days at 2-8 °C.

ASSAY PROCEDURE

Allow the plate to come to room temperature, open it and, if moisture is present on the wells, wait until it has evaporated. Withdraw 5 µL of samples and/or controls and wait until it has been completely absorbed before handling the plate. Cover it tightly and place it in a moist chamber during the incubation time (72 h for IgA and IgG for the procedure 1 or 3, 18 hours for the procedure 2. For IgM 96 hours for the procedure 1 or 3, 24 hours for the procedure 2).

READING TEST RESULTS

Measure the precipitation ring to the nearest 0.1 mm, after the required period according to the procedure should be followed, and the type of protein.

Procedure 1

Reproduce on a plotting paper, the square of the precipitation rings and the concentration of the controls. A line with an intercept should be in the range 10-12 mm². Sample values are determined by interpolation.

Procedure 2

Reproduce on a plotting paper, the square of the precipitation rings against their concentration logarithms. The graph will be a straight line only for a low values. For high values may curve a little. Sample values are obtained by interpolation.

Procedure 3

Read on the enclosed Reference Table, the concentration value corresponding to the precipitation ring diameter. The ring value obtained for the control should have a confidence limit of 0.2 mm (from the reference table).

NOTES

The diffusion time and the reading time depend on the concentration and the specific diffusing protein. After 72 hours or 96 hours, the diffusion of the protein, at any concentration, is completed. For lower concentrations, it is possible to read at lower times (ie 36 h). However, in such cases it is advisable to read again after 3-5 hours. If the diameter is still the same, it is possible to set the concentration, on the contrary, if the diameter is different, the ring should be remeasured after a further 3-5 hours.

QUALITY CONTROL

It is suggested to perform an internal quality control. Not included in the kit.

ANALYTICAL PERFORMANCE

- *Measure's limit.* IgA 70-1050 mg/dL
IgG 300-3500 mg/dL
IgM 40-500 mg/dL

- Precision

IgA mg/dL	Intra-assay		Inter-assay		
	Mean	298.98	453.10	299.62	454.61
SD	5.31	7.07	6.19	6.68	
CV%	1.78	1.56	2.07	1.47	
N	10	10	20	20	
IgG mg/dL	Mean	1299.80	1869.18	1302.32	1863.52
	SD	20.78	26.68	21.90	25.23
SD	1.60	1.43	1.68	1.35	
CV%	10	10	20	20	
IgM mg/dL	Mean	101.00	270.60	102.10	268.60
	SD	2.59	3.37	2.30	4.05
SD	2.57	1.25	2.25	1.51	
CV%	10	10	20	20	

- *Correlation.* This assay (x) was compared with a similar commercial method (y). The results were:

	N = 70	r = 0.97143	y = 1.001x + 2.964
IgA	N = 70	r = 0.99942	y = 0.9835x + 23.161
IgM	N = 70	r = 0.9995	y = 0.989x + 2.17

REFERENCE VALUES

IgA 90-450 mg/dL
IgG 800-1800 mg/dL
IgM 60-280 mg/dL

REFERENCES

Fahey et al., *J. Immunol.*, 94: 84 (1965)
Mancini et al., *Immunochemistry*, 2: 235 (1965)

RI33006-1/0404