STEL 3-Lyser

CONTENTS / CONTENIDO	
REF	3802010 STEL 3-Lyser 1 x 0.5 L
Only for professional <i>in vitro</i> diagnostic use Sólo para uso profesional de diagnóstico <i>in vitro</i>	

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PRINCIPLE

STEL 3-Lyser is used on STEL 3 Auto hematology analyzer for white blood cell counting, differentiating and hemoglobin measurement.

In the process of blood cell counting, red blood cells are dissolved by Lyse and release hemoglobin. Hemoglobin reacts with Lyse; the reaction produces a solution with particular color, so concentration of the solution can be determined by colorimetric method. At the same time, white blood cells release cytoplasm and shrink; size of solid content, such as, nucleus and granules left in different types of white blood cells is different. Based on electrical non-conductivity of blood cells, auto hematology analyzer can classify white blood cells into three groups, such as, lymphocytes, intermediate cells and neutrophils by detecting changes in electrical resistance when suspended particles in electrolyte solution passing through the counting hole.

REAGENT COMPOSITION

Quaternary ammonium salt <10.0%, NaCl <0.15%, Stabilizer <0.12%

STORAGE AND STABILITY

Store at 5-35°C. Relative humidity <90%. Keep the product away from direct sunlight, and do not put it in refrigerator. After opening, the reagent is valid for 30 days. Do not use the diluent if there are any precipitates, particles or floccules in it.

INSTRUCTIONS FOR USE

1. Open the external packaging material, connect pipeline of analyzer with mouth of reagent bottle.
2. When using the product, please refer to instructions for the STEL 3 auto hematology analyzer.

NOTES

1. Dispose the waste and contaminated packing material according to local laws and regulations.
2. Test results may be affected under following conditions:
 - The reagent is decayed or expired;
 - The reagent is contaminated by dust, etc.
3. This reagent is only for in-vitro diagnosis, do not drink it, please put on gloves during operation. In case of contact with skin, eyes and mouth, rinse with plenty of water immediately. See the doctor if necessary.

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