PLASMA PROTEIN MULTICALIBRATOR

CONTENTS / CONTENIDO		
<table border="1"> <tr> <td style="text-align: center;">REF</td> <td>3910005 Plasma Protein Multicalibrator 1 x 2 mL</td> </tr> </table>	REF	3910005 Plasma Protein Multicalibrator 1 x 2 mL
REF	3910005 Plasma Protein Multicalibrator 1 x 2 mL	
<p>Only for <i>in vitro</i> diagnostic use Sólo para uso diagnóstico <i>in vitro</i></p>		

Plasma Protein Multicalibrator Multicalibrador de Proteínas Plasmáticas

Turbidimetry / Turbidimetria

PRINCIPLE

Plasma Protein Multicalibrator is a liquid serum from human origin containing protein concentrations suitable to calibrate determinations of the plasma proteins by turbidimetry method.

REAGENT COMPOSITION

Plasma Protein Multicalibrator.

Human serum. Sodium azide 0.95 g/L.

Precautions. The reagent, from human donors, has been given negative results to anti-HIV 1/2, HBsAg and anti-HCV. Handle cautiously is recommended.

The reagent contains sodium azide. Avoid any contact with skin or mucous.

REAGENT PREPARATION

Plasma Protein Multicalibrator. The Calibrator is ready- to-use.

STORAGE AND STABILITY

- Store at 2-8°C.
The reagent is stable until the expiration date printed on the label. Keep the vial tightly closed at 2-8°C and prevent contaminations during its use.
Do not use the control after the expiration date.
- Presence of particles and turbidity is sign of deterioration.
- Opened vials are stable for 30 days at 2-8°C.

INSTRUCTIONS FOR USE

To be used in turbidimetric assays according to the instructions of each corresponding analytical method.

ASSIGNED VALUES

Concentration values are traceable to the Standard Reference Material CRM 470 (Institute for Reference Materials and Measurements, IRMM).

NOTES

- Bacterial contamination may cause reductions in the stability of many components of the Calibrator.
- When pipetting the Calibrator to the instrument avoid the formation of foam, a common cause of error in the calibration set-up.

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