PLASMA PROTEIN CONTROL N-I

CONTENTS / CONTENIDO	
REF	3915010 Plasma Protein Control N-I 1 x 2 mL
Only for <i>in vitro</i> diagnostic use Sólo para uso diagnóstico <i>in vitro</i>	

Plasma Protein Control N-I Control Proteínas Plasmáticas N-I

Turbidimetry / Turbidimetria

PRINCIPLE

Plasma Protein Control is a human serum containing different concentrations of plasma proteins suitable for the accuracy and precision control of the analytical methods in clinical laboratories.

REAGENT COMPOSITION

Plasma Protein Control N-I.

Human serum. Sodium azide 0.95 g/L.

Precautions. The reagent, from human donors, has been given negative results to anti-HIV 1/2, HBsAg and anti-HCV. Handle cautiously is recommended.

The reagent contains sodium azide. Avoid any contact with skin or mucous.

REAGENT PREPARATION

Plasma Protein Control N-I. The Control is ready-to-use.

STORAGE AND STABILITY

- Store at 2-8°C.
The reagent is stable until the expiration date printed on the label. Keep the control tightly closed at 2-8°C and prevent contaminations during its use.
Do not use the control after the expiration date.
- Presence of particles and turbidity is sign of deterioration.
- Opened vials are stable for 30 days at 2-8°C.

INSTRUCTIONS FOR USE

Use the Control as a patient sample following the instructions of the corresponding analytical method.

ASSIGNED VALUES

The assigned concentration values of each parameter and their intervals have been assigned using the method specified in the table below and are given for orientation only. Results may slightly differ from these values because of the different methods and instruments used in the different laboratories. It is useful that every laboratory establishes its own precision parameters.

NOTES

- Bacterial contamination may cause reductions in the stability of many components of the Control.
- When pipetting the Control to the instrument avoid the formation of foam, a common cause of erroneous results.

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