HbA1c – Calibrator

CONTENTS / CONTENIDO			
REF	3955005	HbA1c Calibrator	4 x 0.25 mL
For <i>in vitro</i> diagnostic use only / Sólo para uso diagnóstico			

HbA1c- Calibrator HbA1c- Calibrador

Latex Turbidimetry / *Turbidimetria Látex*

PRINCIPLE

HbA1c-Calibrator is a group of four liquid and stable calibrators based in human erythrocytes suitable to calibrate determinations of human hemoglobin A1c (HbA1c) by turbidimetry method.

REAGENTS COMPOSITION

CAL 1	Human erythrocytes. Stabilizers.
CAL 2	
CAL 3	
CAL 4	

Precautions. The reagent, from human donors, has been given negative results to anti-HIV 1/2, HBsAg and anti-HCV. It is recommended to handle calibrators with the same precautions as patient samples.

REAGENT PREPARATION

Ready to use.

STORAGE AND STABILITY

1.  The reagents will remain stable until the expiration date printed on the label, when stored tightly closed at 2-8°C and protected from light, heat and contaminations during their use. Do not use the reagents after the expiration date.

PROCEDURE

The calibration curve produced with *HbA1c Calibrator* will be stable for at least 6 days (2-Reagents system) and 8 weeks (3-Reagents system) on most analyzers. The calibrators should be processed with the *HbA1c-Turbidimetric Kit ref: 3155005* in the same way as patient samples. Follow the instructions of the analyzer and Reagents Kit used in the assay for specific instrument calibration process.

LIMITATIONS

Inaccurate results may be resulted by improper pipetting, inadequate mixing, inadequate storing and poorly calibrated analyzers.

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CONCENTRATIONS ACCORDING TO DCCT/NGSP / CONCENTRACIONES CONFORME A LA DCCT/NGSP:

Calibrator / Calibrador	3- Reagents System Sistema 3- Reactivos	2- Reagents System Sistema 2- Reactivos
Level / Nivel 1	5.08%	5.05%
Level / Nivel 2	8.39%	8.45%
Level / Nivel 3	11.5%	12.4%
Level / Nivel 4	14.7%	15.5%

The concentration table is applicable to this lot and sequential sublots, indicated in the label with alphabetic letter following the lot number.

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