# ALBUMIN


Colorimetric method

ENDPOINT

<b>REF CT10020</b> 2 x 40 mL	<b>REF CT10022</b> 6 x 40 mL	<b>REF CT10025</b> 10 x 60 mL
<b>CONTENTS</b> R1.Reagent 2 x 40 mL	<b>CONTENTS</b> R1.Reagent 6 x 40 mL	<b>CONTENTS</b> R1.Reagent 10 x 60 mL
For <i>in vitro</i> diagnostic use only		

## SUMMARY

The method<sup>1</sup> is based on the specific binding of bromocresol green (BCG), an anionic dye, and the protein at acid pH with the resulting shift in the absorption wavelength of the complex. The intensity of the color formed is proportional to the concentration of albumin in the sample.

## REAGENT

### R1. Bromocresol reagent.

Succinate buffer 75 mmol/L pH 4.2, BCG 0.12 mmol/L, tensioactive 2 g/L (w/v).

## PREPARATION

The Reagent is ready-to-use.

## STORAGE AND STABILITY

 Store at 2-8°C. The Reagent is stable until the expiry date stated on the label. After daily use stored tightly closed and protected from light. On board this reagent is stable 30 days. Discard the reagent if the blank presents an absorbance over 0.250 at 630 nm against distilled water.

## SAMPLE COLLECTION

Serum, EDTA plasma unhemolyzed. Albumin in serum and plasma is stable for 2 weeks at 2-8°C, and for up to 4 months at -20°C.

## INTERFERENCES

- Lipemia (inralipid >1,25 g/L) may affect the results.
- Bilirubin (40 mg/dL) does not interfere.
- Hemoglobin (1 g/L) may affect the results.
- Other drugs and substances may interfere<sup>5</sup>.
- Specimens containing dextran should be avoided.

## INSTRUMENTATION AND MATERIALS

- LIDA analyzer.
- Laboratory equipment.
- iso-Clean Solution. Ref. CT18002.
- Multicalibrator CC/H 10x3 mL Ref. CT19750

## AUTOMATED PROCEDURE

A graphic display pictures the specific sets corresponding to the technical application outlined for this test. Any new application, to the instrument should be validated to confirm that results meet the analytical performance of the method. It is recommended to validate periodically the instrument.

## CALIBRATION

Recalibrate weekly, when a new lot of reagent is used, when control recovery falls out of the expected range or when adjustments are made to the instrument. Two point calibration is recommended (S1: NaCl 9 g/L and S2: Calibrator). A reagent blank should be run daily before sample analysis.

## RESULTS

Samples with concentrations higher than 6 g/dL should be diluted 1:2 with saline and assayed again. Multiply the results by 2. If results are to be expressed as SI units apply: g/dL x 10 = g/L

## EXPECTED VALUES<sup>1</sup>

Serum, plasma

Adults	3.81 - 4.65 g/dL (38.1 - 46.5 g/L)
--------	------------------------------------

The range of values for hospitalized individuals varies between 1.4 and 4.8 g/dL.

It is recommended that each laboratory establishes its own reference range.

## QUALITY CONTROL

To ensure adequate quality control (QC), each run should include a set of controls (normal and abnormal) with assayed values handled as unknowns. Each laboratory should establish its own Quality Control scheme and corrective actions if controls do not meet the acceptable tolerances.

## DIAGNOSTIC CHARACTERISTICS<sup>4</sup>

The serum content of the soluble proteins, those circulating in extracellular and intracellular fluids, has been used as a marker to aid in clinical diagnosis. The main diagnostic tests are those measuring serum total protein and serum albumin.

*Hiperproteinemia o hiperalbuminemia*, usually occurs during multiple mieloma caused by high levels of the monoclonal immunoglobulins, dehydration, excessive water loss, as in severe vomiting, diarrhea, Addison's disease or diabetic acidosis. The *hemoconcentration*, decrease in the volume of plasma water, is reflected as a relative hyperproteinemia since concentration of all the individual plasma proteins are increased to the same degree.

*Hypoproteinemia o hypoalbuminemia* usually occurs in edema, malnutrition, nephrotic syndrome, malabsorption and severe liver cirrhosis. Since albumin is present in such high concentration low levels of this protein alone may also cause hypoproteinemia.

Clinical diagnosis should not be made on findings of a single test result, but should integrate both clinical and laboratory data.

## PERFORMANCE CHARACTERISTICS

Performance characteristics are available on request.

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*In vitro diagnostic reagent for the determination of albumin in serum or plasma.*

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## CALIBRATION PARAMETERS

Rule:	Linear (two points)	Calibrator 1:	Deionized Water
Factor K		Calibrator 2:	Multicalibrator
Replicates	1	Calibrator 3:	
Interval		Calibrator 4:	
Sensitivity		Calibrator 5:	
Coefficient		Calibrator 6:	
Error Limit			
Blank response			
Difference Limit			
SD			

## TEST PARAMETERS

No.	
Test	ALB
Method	Endpoint
Direction	Ascend
Unit	g/dL
Decimals	1
Prim. Wave.	630
Sec. Wave.	0
Sample Vol. (µL)	3
R1 Vol. (µL)	400
R2 Vol. (µL)	0
Line. Limit	20
Incubation	0
Reaction	0   5
Antigen	<input type="checkbox"/> Check Antigen
Substrate	0
Response	L Limit 0   U Limit 0
S. Volume	0
Ratio	0
R1 Blank	L Limit 0   U Limit 0
Mix. R Blank	L Limit 0   U Limit 0
Linearity	L Limit 0   U Limit 6
Full Name	ALBUMIN
Print No.	

For other instrument software version please contact with LINEAR