# CALCIUM ARSENAZO III

TOTAL


*Colorimetric method*

ENDPOINT

<b>REF CT10110</b> 2 x 40 mL	<b>REF CT10112</b> 6 x 40 mL
<b>CONTENTS</b> R1.Reagent 2 x 40 mL	<b>CONTENTS</b> R1.Reagent 6 x 40 mL
For <i>in vitro</i> diagnostic use only	

## SUMMARY

The method<sup>1</sup> is based on the specific binding of arsenazo III and calcium at acid pH with the resulting shift in the absorption wavelength of the complex. The intensity of the chromophore formed is proportional to the concentration of total calcium in the sample.

## REAGENT

**R1.** Arsenazo III 120 mmol/L, imidazol 75 mmol/L, pH 6.5.

## PREPARATION

The Reagent is ready-to-use.

## STORAGE AND STABILITY

Store at 2-8°C. The Reagent is stable until the expiry date stated on the label. After daily use stored tightly closed and protected from light. On board this reagent is stable 30 days.

Discard the reagent if the blank presents an absorbance over 0.500 at 650 nm against distilled water.

## SAMPLE COLLECTION

Serum or heparinized plasma and urine (see notes). Other anticoagulants (EDTA, oxalate and citrate) must not be used.

Calcium in serum or plasma is stable for 10 days at 2-8°C. Freeze for longer storage.

Calcium in acidified samples of urine (see Notes) is stable for 10 days at 2-8°C.

## INTERFERENCES

- Lipemia (intralipid >1 g/L) may affect the results.
- Bilirubin (40 mg/dL), hemoglobin (12 g/L), do not interfere.
- Other drugs and substances may interfere<sup>3</sup>.
- Many detergents and water supplies (see Notes).

## INSTRUMENTATION AND MATERIALS

- LIDA analyzer.
- Laboratory equipment.
- iso-Clean Solution. Ref. CT18002.
- Multicalibrator CC/H 10x3 mL Ref. CT19750

## AUTOMATED PROCEDURE

A graphic display pictures the specific sets corresponding to the technical application outlined for this test.

Any new application to the instrument should be validated to confirm that results meet the analytical performance of the method.

It is recommended to validate periodically the instrument.

## CALIBRATION

Recalibrate weekly, when a new lot of reagent is used, when control recovery falls out of the expected range or when adjustments are made to the instrument. Two point calibration is recommended (S1: NaCl 9 g/L and S2: Calibrator). A reagent blank should be run daily before sample analysis.

## RESULTS

Samples with concentrations higher than 20 mg/dL should be diluted 1:2 with saline and assayed again. Multiply the results by 2.

If results are to be expressed as SI units apply: mg/dL x 0.25 = mmol/L.

## EXPECTED VALUES<sup>4</sup>

Serum, plasma

Newborns (< 10 days)	7.6-10.4 mg/dL (1.9-2.6 mmol/L)
Children (2-12 years)	8.8-10.4 mg/dL (2.2-2.6 mmol/L)
Adults (12-60 years)	8.4-10.2 mg/dL (2.1-2.5 mmol/L)
Urine. Adults (normal diet)	100-300 mg/24-h (25-75 mmol/24-h)

It is recommended that each laboratory establishes its own reference range.

## QUALITY CONTROL

To ensure adequate quality control (QC), each run should include a set of controls (normal and abnormal) with assayed values handled as unknowns. Each laboratory should establish its own Quality Control scheme and corrective actions if controls do not meet the acceptable tolerances.

## DIAGNOSTIC CHARACTERISTICS

Calcium exists in the blood in three forms: ionized (13%), complexed (47%) and bound to protein, mainly albumin (40%).

A depressed concentration of total calcium can be due to hypoproteinemia, but the concentration of physiologically active (ionized) calcium in such case may be normal.

Elevated serum calcium levels occur in hyperparathyroidism, vitamin-D poisoning, and sarcoidosis. The plasma level in calcium is greatly affected by the plasma level of inorganic phosphate. In most cases, there is an inverse relationship between calcium and inorganic phosphate.

Conditions associated with *hypercalcemia*, such as primary hyperparathyroidism are usually associated with *hypophosphatemia*; the opposite is true as well. Urine calcium excretion parallels the serum calcium level. Large amounts of calcium are excreted in the urine in hyperparathyroidism, metabolic acidosis, renal tubular insufficiency, multiple myeloma and bone malignancies.

Clinical diagnosis should not be made on findings of a single test result, but should integrate both clinical and laboratory data.

## NOTES

- Most of the detergents and water softening products used in the labs contain chelating agents. A defective rinsing will invalidate the procedure. Keep the glassware acid washed and thoroughly rinsed at all times.
- Collect a 24-hour urine specimen into a plastic bottle containing 10 mL of 50% (v/v) HCl. Centrifuge or filter, and dilute 1:2 with distilled water before testing.

## PERFORMANCE CHARACTERISTICS

Performance characteristics are available on request.

## BIBLIOGRAPHY

1. Michaylova, V. and Illkova, P. Anal. Chim. Acta. 53 : 194 (1971).
2. Friedman and Young. Effects of disease on clinical laboratory tests, 3th ed. AACC Press, 1997.
3. Young DS. Effects of drugs on clinical laboratory tests, 5th ed. AACC Press, 2000.
4. Tietz. N.W. Clinical Guide to Laboratory Tests, 3<sup>rd</sup> Edition. W.B. Saunders Co. Philadelphia, PA. (1995).

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## REFERENCIAS

1. Michaylova, V. y Illkova, P. Anal. Chim. Acta. 53 : 194 (1971).
2. Friedman and Young. Effects of disease on clinical laboratory tests, 3th ed. AACC Press, 1997.
3. Young DS. Effects of drugs on clinical laboratory tests, 5th ed. AACC Press, 2000.
4. Tietz. N.W. Clinical Guide to Laboratory Tests, 3<sup>rd</sup> Edition. W.B. Saunders Co. Philadelphia, PA. (1995).

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## CALIBRATION PARAMETERS

Rule:	Linear (two points)	Calibrator 1:	Deionized Water
Factor K		Calibrator 2:	Multicalibrator
Replicates	1	Calibrator 3:	
Interval		Calibrator 4:	
Sensitivity		Calibrator 5:	
Coefficient		Calibrator 6:	
Error Limit			
Blank response			
Difference Limit			
SD			

## TEST PARAMETERS

No.		
Test	CA	
Method	Endpoint	
Direction	Ascend	
Unit	mg/dL	
Decimals	1	
Prim. Wave.	670	
Sec. Wave.	0	
Sample Vol. (µL)	3	
R1 Vol. (µL)	300	
R2 Vol. (µL)	0	
Line. Limit	20	
Incubation	0	
Reaction	0	15
Antigen	<input type="checkbox"/> Check Antigen	
Substrate	0	
Response	L Limit 0	U Limit 0
S. Volume	0	
Ratio	0	
R1 Blank	L Limit 0	U Limit 0
Mix. R Blank	L Limit	U Limit 0
Linearity	L Limit 0	U Limit 22
Full Name	CALCIUM	
Print No.		

For other instrument software version please contact with LINEAR.