CHLORIDE

Colorimetric method

ENDPOINT

REF CT10130 2 x 40 mL CONTENTS R1, 2 x 40 mL
For <i>in vitro</i> diagnostic use only

SUMMARY

Chloride ions in the sample quantitatively displaces thiocyanate from mercuric thiocyanate. Liberated thiocyanate ion reacts with ferric ion forming a red ferric-thiocyanate complex proportional to the concentration of chloride present in the sample.^{1,2}

REAGENT

R1. Thiocyanate reagent. Mercuric thiocyanate 2 mmol/L, mercuric nitrate 0.1 mmol/L, iron nitrate 30 mmol/L, HNO₃ 45 mmol/L. (see Notes). X_n

PREPARATION Ready-to-use.

STORAGE AND STABILITY

 Store at 2-8°C.

The kit is stable until the expiry date stated on the label. Do not use reagents over the expiration date. Store the vials tightly closed, protected from light and prevented contaminations during the use. On board is stable 30 days.

Discard If appear signs of deterioration:

- Presence of particles and turbidity.
- Blank absorbance (A) at 470 nm > 0.050 in 1cm cuvette.

SAMPLE COLLECTION

Serum, heparinized plasma and CSF. The sera are stable in capped tubes for at least 4 hours at room temperature, 2 days refrigerated (4-8°C) and several months frozen (-20°C). Cerebrospinal fluid should be collected in three sterile tubes, with samples from the first or second tube used for chloride determinations.

INTERFERENCES

- Lipemia (intralipid >1 g/L) may affect the results.
- Bilirubin (40 mg/dL) does not interfere.
- Hemoglobin (12 g/L) does not interfere.
- Other drugs and substances may interfere³.

INSTRUMENTATION AND MATERIALS

- LIDA analyzer.
- Laboratory equipment.
- iso-Clean Solution. Ref. CT18002.
- Multicalibrator CC/H 10x3 mL Ref. CT19750

AUTOMATED PROCEDURE

A graphic display pictures the specific sets corresponding to the technical application outlined for this test.

Any new application, to the instrument should be validated to confirm that results meet the analytical performance of the method.

It is recommended to validate periodically the instrument.

CALIBRATION

Recalibrate weekly, when a new lot of reagent is used, when control recovery falls out of the expected range or when adjustments are made to the instrument. Two point calibration is recommended (S1: NaCl 9 g/L and S2: Calibrator). A reagent blank should be run daily before sample analysis.

RESULTS

Samples with concentrations higher than 125 mEq/L (125 mmol/L) should be diluted 1:2 with distilled water and assayed again. Multiply the results by 2. Concentrations less than 75 mmol/L, should be brought within linear range by increasing the amount of serum used.

EXPECTED VALUES

Serum, plasma	98 - 111 mEq/L (98 - 111 mmol/L)
CSF	120 - 130 mEq/L (120 - 130 mmol/L)

It is recommended that each laboratory establishes its own reference range.

QUALITY CONTROL

To ensure adequate quality control (QC), each run should include a set of controls (normal and abnormal) with assayed values handled as unknowns. Each laboratory should establish its own Quality Control scheme and corrective actions if controls do not meet the acceptable tolerances.

DIAGNOSTIC CHARACTERISTICS

Sodium and chloride represent the majority of the osmotically active constituents of plasma. As a result, chloride is significantly involved in maintenance of water distribution, osmotic pressure, and anion-cation balance in the extracellular fluid compartment.

Hypochloremia (decreased plasma Cl⁻ concentration) is observed in salt-losing nephritis as associated with chronic pyelonephritis. In Addison's disease, Cl⁻ levels as well as Na⁺ levels may drop significantly during the Addisonian crisis and in certain types of metabolic acidosis (e.g., diabetic ketoacidosis and renal failure) and aldosteronism. In metabolic alkalosis, plasma levels of Cl⁻ tend to fall while HCO₃⁻ levels increase.

Hyperchloremia (increased plasma Cl⁻ concentration) occurs with dehydration, renal tubular acidosis, diabetes insipidus, acute renal failure, adrenocortical hyperfunction and metabolic acidosis. Extremely high dietary in take of salt and overtreatment with saline solutions are also causes of hyperchloremia.

Clinical diagnosis should not be made on findings of a single test result, but should integrate both clinical and laboratory data.

PERFORMANCE CHARACTERISTICS

Performance characteristics are available on request.

BIBLIOGRAPHY

1. Schoenfeld R.G., and Lewellen Clin. Chem. 10 : 533 (1964).
2. Tietz. N.W. Clinical Guide to Laboratory Tests, 3rd Edition. W.B. Saunders Co. Philadelphia, PA. (1995).
3. Young DS. Effects of drugs on clinical laboratory tests, 5th ed. AACC Press, 2000.

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REFERENCIAS

1. Schoenfeld R.G.y Lewellen, C.J. Clin.Chem. 10 : 533 (1964).
2. Tietz. N.W. Clinical Guide to Laboratory Tests, 3rd Edition. W.B. Saunders Co. Philadelphia, PA. (1995).
3. Young DS. Effects of drugs on clinical laboratory tests, 5th ed. AACC Press, 2000.

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CALIBRATION PARAMETERS

Rule:	Linear (Endpoint)	Calibrator 1: 0	Deionized Water
Factor K		Calibrator 2:	Chloride
Replicates	3		
Interval			
Sensitivity			
Coefficient			
Error Limit			
Blank response			
Difference Limit			
SD			

TEST PARAMETERS

No.	
Test	CHLORIDE
Method	Endpoint
Direction	Ascend
Unit	mmol/L
Decimals	0
Prim. Wave.	510
Sec. Wave.	0
Sample Vol. (µL)	3
R1 Vol. (µL)	300
R2 Vol. (µL)	
Line. Limit	20
Incubation	0
Reaction	0 25
Antigen	<input type="checkbox"/> Check Antigen
Substrate	0
Response	L Limit 0 U Limit 0
S. Volume	0
Ratio	0
R1 Blank	L Limit 0 U Limit 0
Mix. R Blank	L Limit 0 U Limit 0
Linearity	L Limit 70 U Limit 125
Full Name	CHLORIDE
Print No.	

For other instrument software version please contact with LINEAR