# CREATINE KINASE BR

CK NAC-ACTIVATED




UV enzymatic method

KINETIC

<b>REF CT10160</b> 1 x 50 mL	<b>REF CT10162</b> 3 x 50 mL
<b>CONTENTS</b> R1.Reagent 1 x 40 mL R2.Reagent 1 x 10 mL	<b>CONTENTS</b> R1.Reagent 3 x 40 mL R2.Reagent 3 x 10 mL
For <i>in vitro</i> diagnostic use only	

## SUMMARY

Creatine kinase (CK) catalyzes the reaction between creatine phosphate (CP) and adenosine 5'-diphosphate (ADP) with formation of creatine and adenosine 5'-triphosphate (ATP). The latter phosphorylates glucose to glucose-6-phosphate (G6P) in the presence of hexokinase (HK). G6P is oxidized to Gluconate-6P in the presence of reduced nicotinamide-adenine dinucleotide phosphate (NADP) in a reaction catalyzed by glucose-6-phosphate dehydrogenase (G6P-DH). The conversion is monitored kinetically at 340 nm by the rate of increase in absorbance resulting from the reduction of NADP to NADPH proportional to the activity of CK present in the sample. In this test<sup>1,2</sup> the presence of N-acetylcysteine (NAC) allows the optimal activation of the enzyme.

This test has been formulated according the standardized method described by IFCC. Clin Chem Lab Med 2002; 40(6) : 635-642.

## REAGENTS

**R1. Buffer/Glucose/NAC.** Imidazol buffer 100 mmol/L pH 6.7, glucose 20 mmol/L, NAC 20 mmol/L, magnesium acetate 10 mmol/L, NADP 2.5 mmol/L, HK  $\geq 4$  KU/L, EDTA 2 mmol/L.

**R2. Substrate/Coenzymes.** CP 30 mmol/L, AMP 5 mmol/L, ADP 2 mmol/L, di(adenosine-5') pentaphosphate 10  $\mu$ mol/L, G6P-DH  $\geq 1.5$  KU/L.

## PREPARATION

The Reagents are ready-to-use.

## STORAGE AND STABILITY

 Store at 2-8°C. The Reagents are stable until the expiry date stated on the label. After daily use stored tightly closed and protected from light. On the board the reagents are stable 19 days. Discard the reagent if the blank presents an absorbance over 0.400 at 340 nm against distilled water.

## SAMPLE COLLECTION

Serum. Stable 8 days at 2-8°C or 1 month at -20°C. Chill the samples as rapidly as possible after collection. Moderately or severely hemolyzed specimens are unsatisfactory for testing, as well as plasmas containing EDTA, heparin, citrate, or fluoride since they may produce unpredictable reaction rates.

## INTERFERENCES

- Lipemia (intralipid  $>5$  g/L) may affect the results.
- Bilirubin ( $< 20$  mg/dL), hemoglobin ( $< 10$  g/L), do not interfere.
- Other drugs and substances may interfere<sup>3</sup>.

## INSTRUMENTATION AND MATERIALS

- LIDA analyzer.
- Laboratory equipment.
- iso-Clean Solution. Ref. CT18002.
- Multicalibrator CC/H 10x3 mL Ref. CT19750

## AUTOMATED PROCEDURE

A graphic display pictures the specific sets corresponding to the technical application outlined for this test.

Any new application, to the instrument should be validated to confirm that results meet the analytical performance of the method. It is recommended to validate periodically the instrument.

## CALIBRATION

No calibration is required. Activity is calculated against theoretical factor.

## CALCULATION

Patient values are calculated automatically by theoretical factor.

## RESULTS

Samples with  $\Delta A/\text{min}$  exceeding 0.270 a 340 nm should be diluted 1:10 with saline and assayed again. Multiply the results by 10.

If results are to be expressed as SI units apply: U/L x 16.67 = nkat/L

## EXPECTED VALUES<sup>1</sup>

Serum

Temperature	25°C	30°C	37°C
Men	$\leq 65$ U/L (1083 nkat/L)	$\leq 105$ U/L (1750 nkat/L)	$\leq 174$ U/L (2900 nkat/L)
Women	$\leq 55$ U/L (917 nkat/L)	$\leq 80$ U/L (1334 nkat/L)	$\leq 140$ U/L (2334 nkat/L)
Children	$\leq 94$ U/L (1570 nkat/L)	$\leq 150$ U/L (2500 nkat/L)	$\leq 225$ U/L (3750 nkat/L)

It is recommended that each laboratory establishes its own reference range.

## QUALITY CONTROL

To ensure adequate quality control (QC), each run should include a set of controls (normal and abnormal) with assayed values handled as unknowns. Each laboratory should establish its own Quality Control scheme and corrective actions if controls do not meet the acceptable tolerances.

## DIAGNOSTIC CHARACTERISTICS

Creatine kinase (CK) values are high in patients with myocardial infarction, progressive muscular dystrophy, alcoholic myopathy, and delirium tremens, but normal in patients with hepatitis and other forms of liver disease. The high values in patients with hypothyroidism reflect the muscle changes in this condition. Although CK is found almost exclusively in myocardium, muscle, and brain and early reports suggested it to be an almost specific index of injury of myocardium and muscle, more recent reports indicate that, inexplicably high serum CK values can occur in patients with pulmonary infarction and pulmonary edema. Specificity of CK assay is enhanced measurement of its isoenzymes.<sup>4-6</sup> Clinical diagnosis should not be made on findings of a single test result, but should integrate both clinical and laboratory data.

## PERFORMANCE CHARACTERISTICS

Performance characteristics are available on request.

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*In vitro diagnostic reagent for the determination of creatine kinase (CK) in serum.*

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*Factor included in the application a 37°C*

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QUALITY CONTROL

**Linear Human Multisera Normal CT19800**

**Linear Human Multisera Abnormal CT19850**

### DETERMINACIONES

NUMBER OF TESTS

**200 tests/kit, Code CT10160**

**600 tests/kit, Code CT10162**

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*(dead volume is not taken in consideration)*

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*For additional information read reagent packaging insert.*

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*The reported application is given only as a guideline. Any application to an instrument should be validated to demonstrate that results meet the performance characteristics of the method.*

## CALIBRATION PARAMETERS

Rule:	One-point linear (K Factor)	Calibrator 1:	
Factor K	8095	Calibrator 2:	
Replicates	1	Calibrator 3:	
Interval		Calibrator 4:	
Sensitivity		Calibrator 5:	
Coefficient		Calibrator 6:	
Error Limit			
Blank response			
Difference Limit			
SD			

## TEST PARAMETERS

No.		
Test	CKNAC	
Method	Kinetic	
Direction	Ascend	
Unit	U/L	
Decimals	0	
Prim. Wave.	340	
Sec. Wave.	546	
Sample Vol. (µL)	5	
R1 Vol. (µL)	200	
R2 Vol. (µL)	50	
Line. Limit	20	
Incubation	10	
Reaction	15	30
Antigen	<input type="checkbox"/> Check Antigen	
Substrate	0	
Response	L Limit 0	U Limit 0
S. Volume	0	
Ratio	0	
R1 Blank	L Limit 0	U Limit 0
Mix. R Blank	L Limit 0	U Limit 0
Linearity	L Limit 0	U Limit 900
Full Name	CKNAC	
Print No.		

For other instrument software version please contact with LINEAR.