GGT BR

IFCC

KINETIC

Quantitative determination of γ -GT activity.

REF CT10192 3 x 50 mL	REF CT10195 12 x 50 mL
CONTENTS R1.Reagent 3 x 40 mL R2.Reagent 3 x 10 mL	CONTENTS R1.Reagent 12 x 40 mL R2.Reagent 12 x 10 mL
For professional <i>in vitro</i> diagnostic use only	

SUMMARY

Gamma-glutamyltransferase (γ -GT) catalyzes the transfer of a γ -glutamyl group from γ -glutamyl-3-carboxy-4-nitroanilide to glycylglycine with the formation of L- γ -glutamyl-glycylglycine and 5-amino-2-nitrobenzoate. The amount of 5-amino-2-nitrobenzoate formed is proportional to the enzyme activity present in the sample monitored kinetically at 405 nm.¹

(L- γ -Glutamyl)-glycylglycine + 5-amino-2-nitrobenzoate

REAGENTS

R1. TRIS 133 mmol/L, glycylglycine 138 mmol/L, sodium azide 0.09%.

R2. L- γ -glutamyl-3-carboxy-4-nitroanilide 23 mmol/L, ethylene glycol 90% mmol/L, sodium azide 0.09%.

Warnings and precautions:

Reagent 1. H319 - H335

P264 - 280 - P305 + P351 + P338 - P337 + P313 - P261 -
P271 - P304 + P340 - P312 - P403 + P233 - P405 - P501

Reagent 2. H302 P264 - P270 - P301 + P312 - P330 - P501

Refer to the safety data sheets (SDS) available for download from our website www.linear.es and take the necessary precautions for the use of laboratory reagents.

PREPARATION. The reagents are ready-to-use.

STORAGE AND STABILITY

Store at 2-8°C. Stable until the expiry date stated on the label. After daily use stored tightly closed and protected from light. Do not use reagents over the expiration date. On board stability is 30 days.

Discard if appear signs of deterioration:

- Presence of particles and turbidity.
- Blank absorbance (A) at 405 nm > 1.400 in 1cm cuvette.

SAMPLE COLLECTION

Serum or EDTA plasma free of hemolysis. Fluoride, citrate and oxalate inhibit γ -GT activity². The enzyme in the sample is stable for at least 1 week at 2-8°C and for at least 2 months when frozen.

Only freeze once. Discard contaminated specimens. The samples should be handled as potentially infectious in accordance with Good Manufacturing Practices.

INTERFERENCES⁶

- Lipemia (intralipid >2.5 g/L) may affect the results.
- Bilirubin (> 10 mg/dL) may affect the results.
- Hemoglobin (> 8 g/L) may affect the results.
- Other drugs and substances may interfere⁶.

INSTRUMENTATION AND MATERIALS

- LIDA 500 analyzer.
- Cleaning solutions Ref. 1185010 and 1185015
- Multicalibrator CC/H 8x5 mL Ref. CT19750 (Optional)
- General laboratory equipment.

AUTOMATED PROCEDURE

A graphic display pictures the specific sets corresponding to the technical application outlined for this test.

Any new application, to the instrument should be validated to confirm that results meet the analytical performance of the method.

It is recommended to validate periodically the instrument.

CALIBRATION

With factor. The use of theoretical factor is recommended to calculate the activity.

With calibrator. It is optional, it is only necessary if you do not apply the factor. Calibration must be performed with two points (S1: distilled water and S2: Calibrator). Verify the working reagent blank every day before its use.

It is recommended to calibrate the assay at least every month, following preventive maintenance of the analyser, when using a new reagent lot number, or when quality control results are out of range.

CALCULATION

The analyzer automatically calculates the concentration of each sample by theoretical factor or calibrator.

RESULTS

Samples with $\Delta A/\text{min}$ exceeding 0.200 at 405 nm should be diluted 1:10 with saline and assayed again. Multiply the results by 10.

If results are to be expressed as SI units apply: $U/L \times 16.67 = \text{nkat/L}$

EXPECTED VALUES⁴

Temperature	37°C	30°C	25°C
Men	10-50 U/L (167-834 nKat/L)	7-35 U/L (117-538 nKat/L)	5-25 U/L (83-417 nKat/L)
Women	8-35 U/L (133-583 nKat/L)	6-25 U/L (100-417 nKat/L)	6-25 U/L (83-300 nKat/L)

It is recommended that each laboratory establishes its own reference range.

QUALITY CONTROL

To ensure adequate quality control (QC), each run should include a set of controls (normal and abnormal) Ref. CT19800-CT19850 with assayed values handled as unknowns. Each laboratory should establish its own quality control scheme and corrective actions if controls do not meet the acceptable tolerances.

DIAGNOSTIC CHARACTERISTICS

Gamma-glutamyl transferase is the most sensitive enzymatic indicator available of hepatobiliary disease. The higher increments are found in cases of intrahepatic or posthepatic biliary obstruction, being more sensitive than alkaline phosphatase when detecting obstructive jaundice, colangitis, and cholecystitis. Elevated levels are noted in the sera of patients with alcoholic cirrhosis and from people who are heavy drinkers. The enzyme levels are important in detecting alcohol induced liver disease correlating well with the duration of the drug action.

PERFORMANCE CHARACTERISTICS

Performance characteristics are available on request.

NOTES

1. Clinical diagnosis should not be made on findings of a single test result but should integrate both clinical and laboratory data.
2. The assay procedure must be followed closely, failure to follow the procedure may lead to inaccurate results.
3. To preserve the good functionality of the kit, do not mix different lots neither reagent's remains.
4. Disposal: Used materials should be discarded according to local regulations.

BIBLIOGRAPHY

1. IFCC Primary Reference Procedures for the Measurement of catalytic activity Concentrations of Enzymes at 37°C. part 6. reference procedure for the measurement of catalytic concentration of g-glutamyltransferase. *J. Clin. Chem. Clin. Biochem.* Since 40, (2002).
2. Kryszewski, A. J. J., Neale, G., Whitfield, J. B. B. & Moss, D. W. W. Enzyme changes in experimental biliary obstruction. *Clin. Chim. Acta* 47, 175-182 (1973).
3. Szasz, G., Rosenthal, P. & Fritzsche, W. Die 1-Glutamyl-Transpeptidase-Aktivität im Serum bei hepatobiliären Erkrankungen. *DMW - Dtsch. Medizinische Wochenschrift* 94, 46-61 (1969).
4. Tietz, N.W. Fundamentals of Clinical Chemistry, p.940. W.B. Saunders Co. Philadelphia, PA. (1987).
5. Friedman y Young. Effects of disease on clinical laboratory tests. 5th ed. AACC (Press 2000).
6. Young DS. Effects of drugs on clinical laboratory tests, 5th ed. AACC Press, 2000.

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REFERENCIAS

1. IFCC Primary Reference Procedures for the Measurement of catalytic activity Concentrations of Enzymes at 37°C. part 6. reference procedure for the measurement of catalytic concentratin of g-glutamyltransferase. *J. Clin. Chem. Clin. Biochem./Since* **40**, (2002).
2. Kryszevski, A. J. J., Neale, G., Whitfield, J. B. B. & Moss, D. W. W. Enzyme changes in experimental biliary obstruction. *Clin. Chim. Acta* **47**, 175-182 (1973).
3. Szasz, G., Rosenthal, P. & Fritzsche, W. Die 1-Glutamyl-Transpeptidase-Aktivität im Serum bei hepatobiliären Erkrankungen. *DMW - Dtsch. Medizinische Wochenschrift* **94**, 46-61 (1969).
4. Tietz, N.W. Fundamentals of Clinical Chemistry, p.940. W.B. Saunders Co. Philadelphia, PA. (1987).
5. Friedman y Young. Effects of disease on clinical laboratory tests. 5th ed. AACC (Press 2000).
6. Young DS. Effects of drugs on clinical laboratory tests, 5th ed. AACC Press, 2000.

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