IRON

FERROZINE


Colorimetric method

ENDPOINT

REF CT10262 3 x 50 mL	REF CT10265 12 x 50 mL
CONTENTS R1.Reagent 3 x 40 mL R2.Reagent 3 x 10 mL	CONTENTS R1.Reagent 12 x 40 mL R2.Reagent 12 x 10 mL
For <i>in vitro</i> diagnostic use only	

SUMMARY

The Fe⁺³ bound to serum ferritin once dissociated in a weak-acid medium by Teepol and guanidium chloride, is reduced by hydroxylamine to Fe⁺², forming the ferrous ion a colored complex with FerroZine[®] proportional to the concentration of iron present in the sample.^{1,2}

® Trade mark of Hach Chemical Co. Ames, Iowa.

REAGENTS

R1. Buffer/Reductant. Guanidine chloride 1.0 mol/L, hydroxylamine 0.6 mol/L, acetate buffer 400 mmol/L pH 4.0, Teepol.

R2. Chromogen. FerroZine 40 mmol/L, sodium acetate 400 mmol/L.

PREPARATION

The Reagents are ready-to-use.

STORAGE AND STABILITY

⚠ Store at 2-8°C. The Reagents are stable until the expiry date stated on the label. After daily use stored tightly closed and protected from light. On board the reagents are stable 30 days.

Discard the reagent if the blank presents an absorbance over 0.050 at 560 nm against distilled water.

SAMPLE COLLECTION

Serum or heparinized plasma. Centrifuge specimen as soon as possible after collection. Hemolyzed samples are rejected. Ruptured red cells falsely elevate the serum results.

Iron in serum is stable for 3 weeks at 2-8°C and for about 7 days at 20-25°C. Freeze for longer storage.

INTERFERENCES

- Lipemia (intralipid >1.25 g/L) may affect the results.
- Bilirubin (< 10 mg/dL) does not interfere.
- Hemoglobin may affect the results.
- Other drugs and substances may interfere³.

INSTRUMENTATION AND MATERIALS

- LIDA analyzer.
- Laboratory equipment.
- iso-Clean Solution. Ref. CT18002.
- Multicalibrator CC/H 10x3 mL Ref. CT19750

AUTOMATED PROCEDURE

A graphic display pictures the specific sets corresponding to the technical application outlined for this test.

Any new application, to the instrument should be validated to confirm that results meet the analytical performance of the method.

It is recommended to validate periodically the instrument.

CALIBRATION

Recalibrate weekly, when a new lot of reagent is used, when control recovery falls out of the expected range or when adjustments are made to the instrument. Two point calibration is recommended (S1: NaCl 9 g/L and S2: Calibrator). A reagent blank should be run daily before sample analysis.

RESULTS

Samples with concentrations higher than 1000 µg/dL should be diluted 1:2 with saline and assayed again. Multiply the results by 2.

If results are to be expressed as SI units apply: µg/dL x 0.179 = µmol/L

EXPECTED VALUES⁴

Serum

Men	60-175 µg/dL (10.7-31.3 µmol/L)
Women	50-170 µg/dL (9.0-30.4 µmol/L)

It is recommended that each laboratory establishes its own reference range

QUALITY CONTROL

To ensure adequate quality control (QC), each run should include a set of controls (normal and abnormal) with assayed values handled as unknowns.

Each laboratory should establish its own Quality Control scheme and corrective actions if controls do not meet the acceptable tolerances.

DIAGNOSTIC CHARACTERISTICS

Following intestinal absorption of iron or erythrocyte destruction, iron ions are released into the plasma where they bind to either apotransferrin or apoferritin proteins to form transferrin and ferritin, respectively. The former helps transport iron to bone marrow for *erythropoiesis*; the latter stores iron in tissues, until is needed.

An *increase* in the iron level in plasma due to rapid destruction of erythrocytes or excessive uptake of iron may also lead to iron overload. The latter causes iron deposition disorders in tissue known as *hemosiderosis* or *hemochromatosis*.

Conversely, a *decrease* in the iron level in plasma due to malnutrition or malabsorption may lead to excessive depletion in iron storage, resulting in *anemia* such as iron-deficiency anemia.

Clinical diagnosis should not be made on findings of a single test result, but should integrate both clinical and laboratory data.

PERFORMANCE CHARACTERISTICS

Performance characteristics are available on request.

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In vitro diagnostic reagent for the determination of iron in serum.

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The reported application is given only as a guideline. Any application to an instrument should be validated to demonstrate that results meet the performance characteristics of the method.

CALIBRATION PARAMETERS

Rule:	Linear (two points)	Calibrator 1:	Deionized Water
Factor K		Calibrator 2:	Multicalibrator
Replicates	1	Calibrator 3:	
Interval		Calibrator 4:	
Sensitivity		Calibrator 5:	
Coefficient		Calibrator 6:	
Error Limit			
Blank response			
Difference Limit			
SD			

TEST PARAMETERS

No.	
Test	FE
Method	Endpoint
Direction	Ascend
Unit	µg/dL
Decimals	0
Prim. Wave.	546
Sec. Wave.	700
Sample Vol. (µL)	40
R1 Vol. (µL)	200
R2 Vol. (µL)	50
Line. Limit	20
Incubation	15
Reaction	-1 25
Antigen	<input type="checkbox"/> Check Antigen
Substrate	0
Response	L Limit U Limit
S. Volume	0
Ratio	0
R1 Blank	L Limit 0 U Limit 0
Mix. R Blank	L Limit 0 U Limit 0
Linearity	L Limit 0 U Limit 1000
Full Name	IRON
Print No.	

For other instrument software version please contact with LINEAR.