IRON


CHROMAZUROL
Colorimetric method

ENDPOINT

REF CT10272 3 x 60 mL CONTENTS R1.Reagent 3 x 60 mL CAL. Standard 1 x 3 mL
For <i>in vitro</i> diagnostic use only

PRINCIPLE

The method is based on the properties of Chromazurol S (CAS), a chromogenic iron-binding dye, that under acidic conditions in presence of cetrimide (CTAB) forms an intense purple complex proportional to the concentration of iron present in the sample.

REAGENT

R1. Chromazurol reagent. Acetate buffer 1 mol/L pH 4.7, chromazurol 200 µmol/L, cetrimida 740 µmol/L, Mg²⁺ 2.3 g/L, thiourea 1 mmol/L, Tween 20 0.1 mL/L (v/v).

C R:35/10 S :26-37/39-45

CAL Iron standard. Ferric ion 100 µg/dL (17.9 µmol/L).

REAGENTS PREPARATION

The reagent is ready to use.

STORAGE AND STABILITY

Store at 2-8°C. The reagent is stable until the expiry date stated on the label. On board the reagent is stable 30 days. After daily use store tightly closed and protected from light.

Discard if appear signs of deterioration:

- Presence of particles and turbidity.
- Blank absorbance (A) at 630 nm ≥ 0.400 against distilled water.

SAMPLE COLLECTION

Serum from fasting, early morning samples.

Hemolyzed samples are rejected.

Iron in serum is stable for 3 weeks at 2-8°C. Freeze for longer storage.

INTERFERENCES

- Lipemia (intralipid >1.25 g/L) may affect the results.
- Bilirubin (< 10 mg/dL) does not interfere.
- Hemoglobin may affect the results.
- Other drugs and substances may interfere³.

INSTRUMENTATION AND MATERIALS

- LIDA analyzer.
- Laboratory equipment.
- iso-Clean Solution. Ref. CT18002.

AUTOMATED PROCEDURE

A graphic display pictures the specific sets corresponding to the technical application outlined for this test.

Any new application, to the instrument should be validated to confirm that results meet the analytical performance of the method.

It is recommended to validate periodically the instrument.

CALIBRATION

Use only the calibrator provided with this reagent. The use of calibrators from others sources may not produce accurate results.

Calibration stability on board is kept for 8 days. An abnormal increase of calibration resets is usually indicative of failure in the performance of the reagents. Recalibrate when a new lot of reagents are used, when control recovery falls out range and after any adjustment made to the analyzer.

RESULTS

Samples with concentrations higher than 1000 µg/dL should be diluted 1:2 with saline and assayed again. Multiply the results by 2.

If results are to be expressed as SI units apply: µg/dL x 0.179 = µmol/L

EXPECTED VALUES⁴

Serum, plasma

Men	60-175 µg/dL (10.7-31.3 µmol/L)
Women	50-170 µg/dL (9.0-30.4 µmol/L)

It is recommended that each laboratory establishes its own reference range.

QUALITY CONTROL

To ensure adequate quality control (QC), each run should include a set of controls (normal and abnormal) with assayed values handled as unknowns.

Each laboratory should establish its own Quality Control scheme and corrective actions if controls do not meet the acceptable tolerances.

DIAGNOSTIC CHARACTERISTICS

Following intestinal absorption of iron or erythrocyte destruction, iron ions are released into the plasma where they bind to either apotransferrin or apoferritin proteins to form transferrin and ferritin, respectively. The former helps transport iron to bone marrow for *erythropoiesis*; the latter stores iron in tissues, until is needed.

An *increase* in the iron level in plasma due to rapid destruction of erythrocytes or excessive uptake of iron may also lead to iron overload. The latter causes iron deposition disorders in tissue known as *hemosiderosis* or *hemochromatosis*.

Conversely, a *decrease* in the iron level in plasma due to malnutrition or malabsorption may lead to excessive depletion in iron storage, resulting in *anemia* such as iron-deficiency anemia.

PERFORMANCE CHARACTERISTICS

Performance characteristics are available on request.

BIBLIOGRAPHY

1. Carter, P. Anal. Biochem. 40:450 (1971).
2. Artiss, J.D., Vinogrador, S., and Zak, B. Clin. Biochem. 14:311 (1981).
3. Young DS. Effects of drugs on clinical laboratory tests, 5th ed. AAC Press, 2000.
4. Tietz, N.W., Fundamentals of Clinical Chemistry, p.940. W.B. Saunders Co., Philadelphia, 1987.

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In vitro diagnostic reagent for the determination of iron in serum.

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The reported application is given only as a guideline. Any application to an instrument should be validated to demonstrate that results meet the performance characteristics of the method.

CALIBRATION PARAMETERS

Rule:	Linear (two points)	Calibrator 1:	Deionized Water
Factor K		Calibrator 2:	Iron standard.
Replicates	1	Calibrator 3:	
Interval		Calibrator 4:	
Sensitivity		Calibrator 5:	
Coefficient		Calibrator 6:	
Error Limit			
Blank response			
Difference Limit			
SD			

TEST PARAMETERS

No.		
Test	FECHROM	
Method	Endpoint	
Direction	Ascend	
Unit	µg/dL	
Decimals	0	
Prim. Wave.	630	
Sec. Wave.	0	
Sample Vol. (µL)	7	
R1 Vol. (µL)	250	
R2 Vol. (µL)	0	
Line. Limit	20	
Incubation		
Reaction	0	30
Antigen	<input type="checkbox"/> Check Antigen	
Substrate	0	
Response	L Limit 0	U Limit 0
S. Volume	0	
Ratio	0	
R1 Blank	L Limit 0	U Limit 0
Mix. R Blank	L Limit	U Limit
Linearity	L Limit 0	U Limit 1000
Full Name	IRON CHROMAZUROL	
Print No.		

For other instrument software version please contact with LINEAR.