# TRIGLYCERIDES MR


Enzymatic colorimetric method

ENDPOINT

<b>REF CT10360</b> 2 x 40 mL	<b>REF CT10362</b> 6 x 40 mL	<b>REF CT10365</b> 10 x 60 mL
<b>CONTENTS</b> R1.Reagent 2 x 40 mL	<b>CONTENTS</b> R1.Reagent 6 x 40 mL	<b>CONTENTS</b> R1.Reagent 10 x 60 mL
For <i>in vitro</i> diagnostic use only		

## SUMMARY

The method<sup>1,2</sup> is based on the enzymatic hydrolysis of serum or plasma triglyceride to glycerol and free fatty acids (FFA) by lipoprotein lipase (LPL). The glycerol is phosphorylated by adenosin triphosphate (ATP) in the presence of glycerolkinase (GK) to form glycerol-3-phosphate (G-3-P) and adenosine diphosphate (ADP). G-3-P is oxidized by glycerophosphate oxidase (GPO) to form dihydroxyacetone phosphate (DHAP) and hydrogen peroxide. A red chromogen is produced by the peroxidase (POD) catalyzed coupling of 4-aminoantipyrine (4-AA) and phenol with hydrogen peroxide (H<sub>2</sub>O<sub>2</sub>), proportional to the concentration of triglyceride in the sample.

## REAGENT

**R1.** PIPES buffer 50 mmol/L pH 6.8, LPL ≥ 12 KU/L, GK ≥ 1 KU/L, GPO ≥ 10 KU/L, ATP 2.0 mmol/L, Mg<sup>2+</sup> 40 mmol/L, POD ≥ 2.5 KU/L, 4-AA 0.5 mmol/L, phenol 3 mmol/L, non-ionic tensioactives 2 g/L (w/v). Biocides.

## PREPARATION

The Monoreagent is ready-to-use.

## STORAGE AND STABILITY

Store at 2-8°C. The Reagent is stable until the expiry date stated on the label. After daily use stored tightly closed and protected from light. Discard the reagent if presents an absorbance above 0.150 at 500 nm against distilled water. On board the reagent is stable 30 days.

## SAMPLES

Serum, EDTA or heparinized plasma free of hemolysis. Remove from cells within 2 hours of venipuncture. Analyze samples immediately or refrigerate. Stable for 1 week at 4-8°C.

## INTERFERENCES

- Lipemia (intralipid >2 g/L) may affect the results.
- Bilirubin (20 mg/dL) does not interfere
- Hemoglobin may affect the results.
- Other drugs and substances may interfere<sup>3</sup>.

## INSTRUMENTATION AND MATERIALS

- LIDA analyzer.
- Laboratory equipment.
- iso-Clean Solution. Ref. CT18002.
- Multicalibrator CC/H 10x3 mL. Ref. CT19750

## AUTOMATED PROCEDURE

A graphic display pictures the specific sets corresponding to the technical application outlined for this test.

Any new application, to the instrument should be validated to confirm that results meet the analytical performance of the method.

It is recommended to validate periodically the instrument.

## CALIBRATION

Recalibrate weekly, when a new lot of reagent is used, when control recovery falls out of the expected range or when adjustments are made to the instrument. Two point calibration is recommended (S1: NaCl 9 g/L and S2: Calibrator). A reagent blank should be run daily before sample analysis.

## RESULTS

Samples with concentrations higher than 800 mg/dL should be diluted 1:2 with saline and assayed again. Multiply the results by 2.

If results are to be expressed as SI units apply:

$$\text{mg/dL} \times 0.0113 = \text{mmol/L}$$

## EXPECTED VALUES <sup>4</sup>

Updated clinical values of triglycerides used to classify risk groups.

Triglycerides	Risk Classification
< 150 mg/dL (< 1.70 mmol/L)	Normal
150-199 mg/dL (1.70-2.25 mmol/L)	Borderline/high
200-499 mg/dL (2.26-5.63 mmol/L)	High
≥ 500 mg/dL (≥ 5.65 mmol/L)	Very high

It is recommended that each laboratory establishes its own reference range.

## DIAGNOSTIC CHARACTERISTICS

The plasma level of lipids (triglycerides and cholesterol) and lipid derivatives, especially lipoproteins (HDL and LDL), aids in the diagnosis of many metabolic disorders. An imbalance in the level of lipoproteins in plasma leads to *hyperlipoproteinemia*, a group of disorders that affects lipid and lipoproteins levels in plasma, causing coronary heart disease (CHD) and atherosclerosis. Each type of hyperlipoproteinemia is associated with an abnormal elevation of triglycerides, cholesterol or lipoprotein subfraction.

Prospective studies<sup>4</sup> indicate that elevated triglycerides are also an independent risk for coronary heart disease. The finding that elevated triglycerides are an independent CHD risk factor suggest that some triglyceride-rich lipoproteins are atherogenic. The latter are partially degraded VLDL, commonly called *remnant lipoproteins*. In clinical practice, VLDL cholesterol is the most readily available measure of atherogenic remnant lipoproteins, and as such can be a target of cholesterol-lowering therapy.

Clinical diagnosis should not be made on findings of a single test result, but should integrate both clinical and laboratory data.

## PERFORMANCE CHARACTERISTICS

Performance characteristics are available on request.

## BIBLIOGRAPHY

1. Buccolo G and David, H. Clin. Chem. 19 : 476 (1973).
2. Fossati, R. and Prencipe, L. Clin. Chem. 28: 2077 (1982).
3. Young DS. Effects of drugs on clinical laboratory tests, 5th ed. AACC Press, 2000.
4. SPECIAL REPORT. Executive Summary of the Third Report of the National Cholesterol Education Program (NCEP) Expert Panel on Detection, Evaluation, and Treatment of High Blood Cholesterol in Adults (Adult Treatment Panel III). JAMA. 285 : 2486 (2001).

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## REFERENCIAS

1. Buccolo G y David, H. Clin. Chem. 19 : 476 (1973).
2. Fossati, R. y Prencipe, L. Clin. Chem. 28: 2077 (1982).
3. Young DS. Effects of drugs on clinical laboratory tests, 5th ed. AACC Press, 2000.
4. SPECIAL REPORT. Executive Summary of the Third Report of the National Cholesterol Education Program (NCEP) Expert Panel on Detection, Evaluation, and Treatment of High Blood Cholesterol in Adults (Adult Treatment Panel III). JAMA. 285 : 2486 (2001).

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*In vitro diagnostic reagent for the determination of triglycerides in serum or plasma.*

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**266 tests/kit, Code CT10360**

**800 tests/kit, Code CT10362 2.000 tests/kit, Code CT10365**

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*(dead volume is not taken in consideration)*

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*For additional information read reagent packaging insert.*

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*The reported application is given only as a guideline. Any application to an instrument should be validated to demonstrate that results meet the performance characteristics of the method.*

**CALIBRATION PARAMETERS**

<b>Rule:</b>	<b>Linear (two points)</b>	<b>Calibrator 1:</b>	<b>Deionized Water</b>
<b>Factor K</b>		<b>Calibrator 2:</b>	<b>Multicalibrator</b>
<b>Replicates</b>	<b>1</b>	<b>Calibrator 3:</b>	
<b>Interval</b>		<b>Calibrator 4:</b>	
<b>Sensitivity</b>		<b>Calibrator 5:</b>	
<b>Coefficient</b>		<b>Calibrator 6:</b>	
<b>Error Limit</b>			
<b>Blank response</b>			
<b>Difference Limit</b>			
<b>SD</b>			

**TEST PARAMETERS**

<b>No.</b>		
<b>Test</b>	<b>TRI</b>	
<b>Method</b>	<b>Endpoint</b>	
<b>Direction</b>	<b>Ascend</b>	
<b>Unit</b>	<b>mg/dL</b>	
<b>Decimals</b>	<b>1</b>	
<b>Prim. Wave.</b>	<b>510</b>	
<b>Sec. Wave.</b>	<b>0</b>	
<b>Sample Vol. (µL)</b>	<b>3</b>	
<b>R1 Vol. (µL)</b>	<b>300</b>	
<b>R2 Vol. (µL)</b>	<b>0</b>	
<b>Line. Limit</b>	<b>20</b>	
<b>Incubation</b>	<b>0</b>	
<b>Reaction</b>	<b>0</b>	<b>25</b>
<b>Antigen</b>	<input type="checkbox"/> <b>Check Antigen</b>	
<b>Substrate</b>	<b>0</b>	
<b>Response</b>	<b>L Limit 0</b>	<b>U Limit 0</b>
<b>S. Volume</b>	<b>0</b>	
<b>Ratio</b>	<b>0</b>	
<b>R1 Blank</b>	<b>L Limit 0</b>	<b>U Limit 0</b>
<b>Mix. R Blank</b>	<b>L Limit 0</b>	<b>U Limit 0</b>
<b>Linearity</b>	<b>L Limit 0</b>	<b>U Limit 800</b>
<b>Full Name</b>	<b>TRIGLYCERIDES</b>	
<b>Print No.</b>		

For other instrument software version please contact with LINEAR