β_2 -microglobulin Turbidimetric

Latex Turbidimetry

REF CT31160 1 x 50 mL CONTENTS R1. Reagent 1 x 40 mL R2. Reagent 1 x 10 mL	REF KR31162 3 x 50 mL CONTENTS R1. Reagent 3 x 40 mL R2. Reagent 3 x 10 mL
For <i>in vitro</i> diagnostic use only	

SUMMARY

The latex particles coated with anti-human β_2 -microglobulin are agglutinated when they react with samples that containing β_2 -microglobulin (β_2 -M). The latex particles agglutination is proportional to the concentration of the β_2 -M in the sample and can be measured by turbidimetry¹.

REAGENTS

R1 Diluent. Tris buffer, 20 mmol/L, pH 8.2.

R2. Latex. Latex particles coated with goat IgG anti-human β_2 -M, pH 8.2.

Optional: β_2 -M Calibrator Set 2 x 1 mL (ref. CT39160).

CAL-S Serum Calibrator. Human serum.

CAL-U Urine Calibrator. Human urine.

For additional information, read β_2 -M Calibrator Set.

Precautions: The reagents contain sodium azide 0.95 g/L. Avoid any contact with skin or mucous.

PREPARATION

R1. Ready to use.

R2. Ready to use. Mix gently the vial by inversion before use (Note 4).

STORAGE AND STABILITY

-  The reagents will remain stable until the expiration date printed on the label. For optimal stability store **tightly closed** at 2-8°C. Do not use the reagents after the expiration date.
- Reagent deterioration: Presence of particles and turbidity.

SAMPLE COLLECTION

Fresh serum. Stable 7 days at 2-8°C or 3 months at -20°C. Samples with presence of fibrin should be centrifuged before testing. Hemolyzed or contaminated samples are not suitable for testing.

Fresh urine. Adjust urine pH at 7.6 by the addition of 0.4 mol/L K_2HPO_4 .

Stable 2 days at 2-8°C or 2 months at -20°C.

INTERFERENCES

- Serum:** Bilirubin (40 mg/dL), hemoglobin (12 g/L), lipemia (> 20 g/L) and rheumatoid factors (300 IU/mL), do not interfere.
 - Urine:** Hemoglobin (12 g/L), creatinine (300 mg/L), uric acid (500 mg/L) and urea (100 mg/L), do not interfere.
- Other substances may interfere⁹.

ADDITIONAL EQUIPMENT

- LIDA Analyzer.
- Laboratory equipment.
- Iso-Clean Solution Ref. CT18002.
- Linear β_2 -M Calibrator Set 2 x 1 mL (ref. CT39160).

AUTHOMATIC TECHNIQUE

For automatic assays the application method is included. Any application to an instrument should be validated to demonstrate that results meet the performance characteristics of the method. It is recommended to validate periodically the instrument.

CALIBRATION

Recalibrate weekly, when a new lot of reagent is used, when control recovery falls out of the expected range or when adjustments are made to the instrument. A reagent blank should be run daily before sample analysis.

EXPECTED VALUES⁷

Serum⁸: 1 – 3 mg/L.

Urine⁸: < 0.3 mg/L.

Each laboratory should establish its own reference range.

QUALITY CONTROL

Control sera are recommended to monitor the performance of manual and automated assay procedures. It is recommended to use Plasma Protein Control N-I (ref: 3915010) and N-II (ref: 3915015).

Each laboratory should establish its own Quality Control scheme and corrective actions if controls do not meet the acceptable tolerances.

DIAGNOSTIC CHARACTERISTICS⁴⁻⁶

β_2 -M is a low molecular weight polipeptide (11,800) that is a constituent of HLA⁴. Because approximately 50% of β_2 -M is produced by lymphocytes and is freely filtered through glomerular basement membrane, β_2 -M levels rise in serum in the presence of glomerular impairment, or lymphocyte activation. Therefore, elevated levels of β_2 -M in serum have been reported in disorders such as chronic inflammatory disorders (rheumatoid arthritis, systemic lupus erythematosus, sarcoidosis), viral infection (Epstein Barr virus, cytomegalovirus, influenza A, HIV), lymphoid malignancies (leukemia, malignant lymphoma, myeloma) and cancer (breast, lung, stomach)⁴⁻⁶.

As for β_2 -M levels in urine, only a trace concentration is detected under normal conditions because up to 99% of β_2 -M passed through glomerular is reabsorbed at renal proximal tubules. Therefore β_2 m levels rise significantly with renal tubular dysfunction, or when serum levels are above. The tubular reabsorption threshold (usually > 4 mg/mL)⁷.

PERFORMANCE CHARACTERISTICS

Performance characteristics are available on request.

NOTES

- This method may be used with different instruments. Any application to an instrument should be validated to demonstrate that results meet the performance characteristics of the method. It is recommended to validate periodically the instrument. Contact to the distributor for any question on the application method.
- The linearity limit depends on the sample/reagent ratio, as well as the analyzer used. It will be higher by decreasing the sample volume, although the sensitivity of the test will be proportionally decreased.
- Clinical diagnosis should not be made on findings of a single test result, but should integrate both clinical and laboratory data.
- For automatic instruments, avoid the presence of bubbles in the reagents that may interfere with the assay results.

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