RF - Turbidimetric

Latex Turbidimetry

| REF CT31300 1 x 50 mL CONTENTS | REF CT31312 3 x 50 mL CONTENTS | REF CT31315 12 x 50 mL CONTENTS |
|--|---|--|
| R1.Reagent 1x 40 mL R2. Reagent 1 x 10 mL | R1.Reagent 3 x 40 mL R2. Reagent 3 x 10 mL | R1.Reagent 12 x 40 mL R2.Reagent 12 x 10 mL |
| For <i>in vitro</i> diagnostic use only | | |

SUMMARY

The latex particles coated with human gammaglobulin are agglutinated when they react with samples that contain rheumatoid factors (RF). The latex particles agglutination is proportional to the concentration of the RF in the sample and can be measured by turbidimetry.^{8,9}

REAGENTS

R1. Diluent. Tris buffer, 20 mmol/L, pH 8.2.

R2. Latex. Latex particles coated with human gammaglobulin, pH 8.2.

Precautions: The reagents contain sodium azide 0.95 g/L. Avoid any contact with skin or mucous.

PREPARATION

R1. Ready to use.

R2. Ready to use. Shake gently the vial before use.

Calibration curve: Prepare two fold dilutions of the Calibrator (1/1, 1/2, 1/4, 1/8) using NaCl 9 g/L as diluent. Use NaCl 9 g/L as Calibrator 0. Calculate the concentration of each RF Calibrator dividing the Calibrator concentration by the corresponding dilution factor.

STORAGE AND STABILITY

-  The reagents will remain stable until the expiration date printed on the label. For optimal stability store **tightly closed** at 2-8°C. Do not use the reagents after the expiration date.
- Reagent deterioration: Presence of particles and turbidity.

SAMPLE COLLECTION

Fresh serum. Stable 7 days at 2-8°C or 3 months at -20°C. Samples with presence of fibrin should be centrifuged before testing. Hemolyzed or contaminated samples are not suitable for testing.

INTERFERENCES

- Bilirubin (40 mg/dL), hemoglobin (4 g/L), lipemia (5 g/L), do not interfere. Other substances may interfere¹¹.

ADDITIONAL EQUIPMENT

- LIDA Analyzer.
- Laboratory equipment.
- iso-Clean Solution. Ref. CT18002.
- Linear RF Calibrator Ref. 3931305.

AUTHOMATIC TECHNIQUE

For automatic assays the application method is included. Any application to an instrument should be validated to demonstrate that results meet the performance characteristics of the method. It is recommended to validate periodically the instrument.

CALIBRATION

Recalibrate weekly, when a new lot of reagent is used, when control recovery falls out of the expected range or when adjustments are made to the instrument. A reagent blank should be run daily before sample analysis.

EXPECTED VALUES¹⁰

Adults: Up to 30 IU/mL.
Each laboratory should establish its own reference range.

QUALITY CONTROL

To ensure adequate quality control (QC), each run should include a set of controls (normal and abnormal) with assayed values handled as unknown. Each laboratory should establish its own Quality Control scheme and corrective actions if controls do not meet the acceptable tolerances.

DIAGNOSTIC CHARACTERISTICS¹⁻⁷

Rheumatoid factors are a group of IgM antibodies directed to determinants in the Fc portion of the immunoglobulin G molecule. Although rheumatoid factors are found in a number of rheumatoid disorders, such as systemic lupus erythematosus (SLE) and Sjögren's syndrome, as well as in nonrheumatic conditions, its central role in clinic lies its utility as an aid in the diagnosis of rheumatoid arthritis (RA).

The rheumatoid factors are found in 70-100% of cases of definite rheumatoid arthritis depending on the test procedure used to detect them.

PERFORMANCE CHARACTERISTICS

Performance characteristics are available on request.

NOTES

- The linearity limit depends on the sample/reagent ratio, as well as the analyzer used. It will be higher by decreasing the sample volume, although the sensitivity of the test will be proportionally decreased.
- Clinical diagnosis should not be made on findings of a single test result, but should integrate both clinical and laboratory data.

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