α_1 -ac. glycoprotein at

Turbidimetric method

REF CT31350 2 x 40 mL CONTENTS R1.Reagent 2 x 40 mL	REF CT31352 6 x 40 mL CONTENTS R1.Reagent 6 x 40 mL
For <i>in vitro</i> diagnostic use only	

SUMMARY

α_1 -acid glycoprotein at is a quantitative turbidimetric assay¹ for the measurement of α_1 -acid glycoprotein (α_1 -AG) in human serum or plasma. Anti-human α_1 -AG antibodies form insoluble complexes when mixed with samples containing α_1 -AG. The light scattering of the immunocomplexes depends upon the α_1 -AG concentration of the patient sample, and can be quantified by comparison from a calibrator of known α_1 -AG concentration.

REAGENTS

R1 α_1 -acid glycoprotein at. Goat antibodies anti-human α_1 -acid glycoprotein, tris buffer 20 mmol/L, pH 8.2. Sodium azide 0.95 g/L.

Precautions: The reagent contains sodium azide 0.95 g/L. Avoid any contact with skin or mucous.

PREPARATION

R1. Ready to use.

Calibration curve:

Prepare two fold dilutions of the Calibrator (1/1, 1/2, 1/4, 1/8) using NaCl 9 g/L as diluent. Use NaCl 9 g/L as Calibrator 0. Calculate the concentration of each α_1 -acid glycoprotein Calibrator dividing the Calibrator concentration by the corresponding dilution factor.

STORAGE AND STABILITY

1.  Store at 2-8°C.

The reagent is stable until the expiration date stated on the label, when stored **tightly closed** at 2-8°C and contaminations are prevented during their use. Do not use the reagent after the expiration date.

2. Presence of particles, turbidity and the absorbance of blank reagent > 0.3 at 340 nm are signs of deterioration.

SAMPLE COLLECTION

Fresh serum or, EDTA or heparinized plasma. Stable 7 days at 2-8°C or 3 months at -20°C.

α_1 -AG in serum or plasma is stable 7 days at 2-8°C or 3 months at -20°C. Samples with presence of fibrin should be centrifuged before testing. Hemolyzed or contaminated samples are not suitable for testing.

INTERFERENCES

- Bilirubin (10 mg/dL), hemoglobin (2 g/L), rheumatoid factors (200 IU/mL), do not interfere.
- Lipemia (≥ 1.25 g/L), may affect the results. Other substances may interfere^{6,7}.

ADDITIONAL EQUIPMENT

- Analyzer LIDA.

AUTHOMATIC TECHNIQUE

For automatic assays the application method is included. Any application to an instrument should be validated to demonstrate that results meet the performance characteristics of the method. It is recommended to validate periodically the instrument.

CALIBRATION

Recalibrate weekly, when a new lot of reagent is used, when control recovery falls out of the expected range or when adjustments are made to the instrument.

A reagent blank should be run daily before sample analysis.

REFERENCE VALUES

Adults²: 50 – 120 mg/dL.

It is recommended that each laboratory establishes its own reference range.

QUALITY CONTROL

To ensure adequate quality control (QC), each run should include a set of controls (normal and abnormal) with assayed values handled as unknown. Each laboratory should establish its own Quality Control scheme and corrective actions if controls do not meet the acceptable tolerances.

DIAGNOSTIC CHARACTERISTICS³⁻⁵

α_1 - acid glycoprotein (orosomucoide) is an acute phase protein³ shown to exert a regulatory, dampening influence on the inflammatory cascade, thereby protecting against tissue damage from excessive inflammation.

α_1 - acid glycoprotein is synthesized by hepatic parenchymal cells, but granulocytes and monocytes also contribute significantly to plasma levels in sepsis. It is classified as a lipocalin, a group of proteins that binds lipophilic substances as progesterone and related hormones⁵.

Plasma concentration increases 3 to 4-fold in most conditions associated with inflammation or tissue necrosis, and may be one of the most reliable indicators of clinical activity of ulcerative colitis⁴. Levels also are increased by glucocorticoid effect, either endogenous (Cushing's syndrome) or exogenous, along with heptaglobin and prealbumin levels. Synthesis and plasma levels are decreased by estrogens⁵.

Patients with proteinuria, α_1 -acid glycoprotein is preferentially excreted in the urine.

PERFORMANCE CHARACTERISTICS

Performance characteristics are available on request.

NOTES

1. The linearity limit depends on the sample/reagent ratio, as well as the analyzer used. It will be higher by decreasing the sample volume, although the sensitivity of the test will be proportionally decreased.
2. Clinical diagnosis should not be made on findings of a single test result, but should integrate both clinical and laboratory data.

BIBLIOGRAPHY

1. Price CP et al. *Ann Clin Biochem* 20: 1-14 (1983).
2. Dati F et al. *Eur J Clin Chem Clin Biochem* 34:517-520 (1996).
3. Kushner, Mackiewicz A. *Disease Markers* 1-11 (1987).
4. Birger K et al. *Scand J Gastroent* 11: 177-183 (1976).
5. Tietz *Textbook of Clinical Chemistry*, 3rd Ed. Burtis CA, Ashwood ER. WB Saunders Co., (1999).
6. Young DS. *Effects of drugs on clinical laboratory tests*. 3th ed. AACC Press (1997).
7. Friedman and Young. *Effects of the disease on clinical laboratory tests*, 3th ed. AACC Press, 1997.

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REFERENCIAS

1. Price CP et al. Ann Clin Biochem 20: 1-14 (1983).
2. Dati F et al. Eur J Clin Chem Clin Biochem 34:517-520 (1996).
3. Kushner, Mackiewicz A. Disease Markers 1-11 (1987).
4. Birger K et al. Scand J Gastroent 11: 177-183 (1976).
5. Tietz Textbook of Clinical Chemistry, 3rd Ed. Burtis CA, Ashwood ER. WB Saunders Co., (1999).
6. Young DS. Effects of drugs on clinical laboratory tests. 3th ed. AACC Press (1997).
7. Friedman and Young. Effects of the disease on clinical laboratory tests, 3th ed. AACC Press, 1997.

Aplicaciones en LIDA 300
Applications on LIDA 300

α_1 - Acid Glycoprotein at

α_1 - Acid Glycoprotein at (2 x 40 mL), Code CT31350
 α_1 - Acid Glycoprotein at (6 x 40 mL), Code CT31352

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CALIBRATION PARAMETERS

Rule:	Logistic 5 P	Calibrator 1:	Deionized Water
Factor K		Calibrator 2:	*PPM1
Replicates	2	Calibrator 3:	*PPM2
Interval		Calibrator 4:	*PPM3
Sensitivity		Calibrator 5:	*PPM4
Error Limit		*Plasma Protein Multicalibrator	
Blank response			
Difference Limit			
SD			

TEST PARAMETERS

No.		
Test	GLYAc	
Method	Endpoint	
Direction	Ascend	
Unit	mg/dL	
Decimals	1	
Prim. Wave.	340	
Sec. Wave.	0	
Sample Vol. (μL)	3	
R1 Vol. (μL)	450	
R2 Vol. (μL)	0	
Line. Limit	20	
Incubation	0	
Reaction	0	20
Antigen	<input type="checkbox"/> Check Antigen	
Substrate	0	
Response	L Limit 0	U Limit 0
S. Volume	0	
Ratio	0	
R1 Blank	L Limit 0	U Limit 0
Mix. R Blank	L Limit 0	U Limit 0
Linearity	L Limit 2	U Limit 250
Full Name	Glycoprotein Acid	
Print No.		

For other instrument software version please contact with LINEAR