C3 at

Turbidimetric method

REF CT31400 2 x 40 mL CONTENTS R1.Reagent 2 x 40 mL	REF CT31402 6 x 40 mL CONTENTS R1.Reagent 6 x 40 mL
For <i>in vitro</i> diagnostic use only	

SUMMARY

C3 at is a quantitative turbidimetric assay^{1,2} for the measurement of (C3) in human serum or plasma.

Anti-human C3 antibodies form insoluble complexes when mixed with samples containing C3. The scattering light of the immunocomplexes depends of the C3 concentration in the patient sample, and can be quantified by comparison from a calibrator of known C3 concentration.

REAGENTS

R1 C3 at. Goat antibodies anti-human C3, tris buffer 20 mmol/L, pH 8.2.

Precautions: The reagent contains sodium azide 0.95 g/L. Avoid any contact with skin or mucous.

PREPARATION

R1. Ready to use.

Calibration curve:

Prepare two fold dilutions of the Calibrator (1/1, 1/2, 1/4, 1/8) using NaCl 9 g/L as diluent. Use NaCl 9 g/L as Calibrator 0. Calculate the concentration of each C3 Calibrator dividing the Calibrator concentration by the corresponding dilution factor.

STORAGE AND STABILITY

1.  The reagents will remain stable until the expiration date printed on the label. For optimal stability store **tightly closed** at 2-8°C. Do not use the reagents after the expiration date.
2. Presence of particles, turbidity and/or the absorbance of blank reagent > 0.3 at 340 nm are sign of deterioration.

SAMPLE COLLECTION

Fresh serum and EDTA or heparinized plasma.

C3 in serum or plasma is stable 7 days at 2-8°C or 3 months at -20°C. Samples with presence of fibrin should be centrifuged before testing. Highly hemolyzed or lipemic samples are not suitable for testing.

INTERFERENCES

- Bilirubin (10 mg/dL), hemoglobin (4 g/L), and rheumatoid factors (300 IU/mL) may affect the results. Lipemia (10 g/L) does not interfere. Other substances may interfere⁵.

ADDITIONAL EQUIPMENT

- Analyzer LIDA.

AUTHOMATIC TECHNIQUE

For automatic assays the application method is included.

Any application to an instrument should be validated to demonstrate that results meet the performance characteristics of the method.

It is recommended to validate periodically the instrument.

CALIBRATION

Recalibrate weekly, when a new lot of reagent is used, when control recovery falls out of the expected range or when adjustments are made to the instrument. A reagent blank should be run daily before sample analysis.

REFERENCE VALUES

Adults³: 90 – 180 mg/dL

Newborn⁴: 70 – 196 mg/dL

It is recommended that each laboratory establishes its own reference range.

QUALITY CONTROL

To ensure adequate quality control (QC), each run should include a set of controls (normal and abnormal) with assayed values handled as unknown. Each laboratory should establish its own Quality Control scheme and corrective actions if controls do not meet the acceptable tolerances.

DIAGNOSTIC CHARACTERISTICS³⁻⁵

C3 is the complement component of highest concentration in plasma and the most functionally linked between the classical and alternative pathway activations¹. Hepatic cells synthesize C3, although bacterial endotoxins induce synthesis through monocytes and fibroblasts.

Increased and decreased levels of C3 both have clinical significance.

Increased levels are closely related with acute-phase response (trauma, inflammatory process), biliary obstruction and focal glomeruloesclerosis².

Decreased levels are related with genetic deficiency (risk for infection, particularly by encapsulated bacteria), or acquired deficiency (collagen vascular diseases and severe infections)⁴.

PERFORMANCE CHARACTERISTICS

Performance characteristics are available on request.

NOTES

1. The linearity limit depends on the sample/reagent ratio, as well as the analyzer used. It will be higher by decreasing the sample volume, although the sensitivity of the test will be proportionally decreased.
2. Clinical diagnosis should not be made on findings of a single test result, but should integrate both clinical and laboratory data.

BIBLIOGRAPHY

1. Lambris John D. *Immunology Today* 9: 387-393 (1988).
2. Price CP et al. *Ann Clin Biochem* 20: 1-14 (1983).
3. Dati F et al. *Eur J Clin Chem Clin Biochem* 34 : 517-520 (1996).
4. Tietz *Textbook of Clinical Chemistry*, 3rd Ed. Burtis CA, Ashwood ER. WB Saunders Co., (1999).
5. Young DS. *Effects of drugs on clinical laboratory tests*. 3th ed. AACC Press (1997).
6. Friedman and Young. *Effects of the disease on clinical laboratory tests*, 3th ed. AACC Press, 1997.

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REFERENCIAS

1. Lambris John D. *Immunology Today* 9: 387-393 (1988).
2. Price CP et al. *Ann Clin Biochem* 20: 1-14 (1983).
3. Dati F et al. *Eur J Clin Chem Clin Biochem* 34:517-520 (1996).
4. Tietz *Textbook of Clinical Chemistry*, 3rd Ed. Burtis CA, Ashwood ER. WB Saunders Co., (1999).
5. Young DS. *Effects of drugs on clinical laboratory tests*. 3th ed. AACC Press (1997).
6. Friedman and Young. *Effects of the disease on clinical laboratory tests*, 3th ed. AACC Press, 1997.

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CALIBRATION PARAMETERS

Rule:	Logistic 5 P	Calibrator 1:	Deionized Water
Factor K		Calibrator 2:	*PPM1
Replicates	2	Calibrator 3:	*PPM2
Interval		Calibrator 4:	*PPM3
Sensitivity		Calibrator 5:	*PPM4
Error Limit		*Plasma Protein Multicalibrator	
Blank response			
Difference Limit			
SD			

TEST PARAMETERS

No.	
Test	C3
Method	Endpoint
Direction	Ascend
Unit	mg/dL
Decimals	1
Prim. Wave.	340
Sec. Wave.	0
Sample Vol. (µL)	3
R1 Vol. (µL)	300
R2 Vol. (µL)	0
Line. Limit	
Incubation	0
Reaction	0 10
Antigen	<input type="checkbox"/> Check Antigen
Substrate	0
Response	L Limit 0 U Limit 0
S. Volume	0
Ratio	0
R1 Blank	L Limit 0 U Limit 0
Mix. R Blank	L Limit 0 U Limit 0
Linearity	L Limit 4 U Limit 500
Full Name	C3
Print No.	

For other instrument software version please contact with LINEAR