C4 at

Turbidimetric method

REF CT31500 2 x 40 mL	REF CT31502 6 x 40 mL
CONTENTS R1.Reagent 2 x 40 mL	CONTENTS R1.Reagent 6 x 40 mL
For <i>in vitro</i> diagnostic use only	

SUMMARY

C4 at is a quantitative turbidimetric assay^{1,2} for the measurement of (C4) in human serum or plasma.

Anti-human C4 antibodies form insoluble complexes when mixed with samples containing C4. The scattering light of the immunocomplexes depends of the C4 concentration in the patient sample, and can be quantified by comparison from a calibrator of known C4 concentration.

REAGENTS

R1 C4 at. Goat antibodies anti-human C4, tris buffer 20 mmol/L, pH 8.2.

Precautions: The reagent contains sodium azide 0.95 g/L. Avoid any contact with skin or mucous.

PREPARATION

R1. Ready to use.

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Prepare two fold dilutions of the Calibrator (1/1, 1/2, 1/4, 1/8) using NaCl 9 g/L as diluent. Use NaCl 9 g/L as Calibrator 0. Calculate the concentration of each C4 Calibrator dividing the Calibrator concentration by the corresponding dilution factor.

STORAGE AND STABILITY

1.  The reagents will remain stable until the expiration date printed on the label. For optimal stability store **tightly closed** at 2-8°C. Do not use the reagents after the expiration date.
2. Presence of particles, turbidity and/or the absorbance of blank reagent > 0.3 at 340 nm are sign of deterioration.

SAMPLE COLLECTION

Fresh serum and EDTA or heparinized plasma.

C4 in serum or plasma is stable 7 days at 2-8°C or 3 months at -20°C. Samples with presence of fibrin should be centrifuged before testing. Highly hemolyzed or lipemic samples are not suitable for testing.

INTERFERENCES

- Bilirubin (10 mg/dL) and rheumatoid factors (400 IU/mL) do not interfere. Hemoglobin (4 g/L) and lipemia (2.5 g/L) may affect the results. Other substances may interfere³.

ADDITIONAL EQUIPMENT

- Analyzer LIDA.

AUTHOMATIC TECHNIQUE

For automatic assays the application method is included.

Any application to an instrument should be validated to demonstrate that results meet the performance characteristics of the method.

It is recommended to validate periodically the instrument.

CALIBRATION

Recalibrate weekly, when a new lot of reagent is used, when control recovery falls out of the expected range or when adjustments are made to the instrument. A reagent blank should be run daily before sample analysis.

REFERENCE VALUES

Adults³: 10 – 40 mg/dL

Newborn⁴: 13 – 38 mg/dL

It is recommended that each laboratory establishes its own reference range.

QUALITY CONTROL

To ensure adequate quality control (QC), each run should include a set of controls (normal and abnormal) with assayed values handled as unknown. Each laboratory should establish its own Quality Control scheme and corrective actions if controls do not meet the acceptable tolerances.

DIAGNOSTIC CHARACTERISTICS³⁻⁵

C4 is the complement component essential for classical pathway activation. Most individuals with C4 deficiency do not have problems with infection, suggesting that the alternative pathway can compensate for the lack in the classical pathway activation in removal of bacterial agents. Hepatic cells synthesize C4, although in less proportion may be synthesized by monocytes and other tissues.

Increased and decreased levels of C4 both have clinical significance.

Increased levels are closely related with acute-phase response (trauma, inflammatory process).

Decreased levels are related with genetic deficiency (autoimmune or collagen vascular disease, particularly Systemic Lupus Erythematosus), or acquired deficiency as a consequence of the consumption in immunocomplexes formation, autoimmune hemolytic anemia and autoimmune nephritis.

PERFORMANCE CHARACTERISTICS

Performance characteristics are available on request.

NOTES

1. The linearity limit depends on the sample/reagent ratio, as well as the analyzer used. It will be higher by decreasing the sample volume, although the sensitivity of the test will be proportionally decreased.
2. Clinical diagnosis should not be made on findings of a single test result, but should integrate both clinical and laboratory data.

BIBLIOGRAPHY

1. Sahu A, Lambris JD. *Immunological Reviews* 180: 35-48 (2001).
2. Price CP et al. *Ann Clin Biochem* 20: 1-14 (1983).
3. Dati F et al. *Eur J Clin Chem Clin Biochem* 34 : 517-520 (1996).
4. Tietz *Textbook of Clinical Chemistry*, 3rd Ed. Burtis CA, Ashwood ER. WB Saunders Co., (1999).
5. Young DS. *Effects of drugs on clinical laboratory tests*. 3th ed. AACC Press (1997).
6. Friedman and Young. *Effects of the disease on clinical laboratory tests*, 3th ed. AACC Press, 1997.

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REFERENCIAS

1. Sahu A, Lambris JD. *Immunological Reviews* 180: 35-48 (2001).
2. Price CP et al. *Ann Clin Biochem* 20: 1-14 (1983).
3. Dati F et al. *Eur J Clin Chem Clin Biochem* 34:517-520 (1996).
4. Tietz *Textbook of Clinical Chemistry*, 3rd Ed. Burtis CA, Ashwood ER. WB Saunders Co., (1999).
5. Young DS. *Effects of drugs on clinical laboratory tests*. 3th ed. AACC Press (1997).
6. Friedman and Young. *Effects of the disease on clinical laboratory tests*, 3th ed. AACC Press, 1997.

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CALIBRATION PARAMETERS

Rule:	Logistic 5 P	Calibrator 1:	Deionized Water
Factor K		Calibrator 2:	*PPM1
Replicates	2	Calibrator 3:	*PPM2
Interval		Calibrator 4:	*PPM3
Sensitivity		Calibrator 5:	*PPM4
Error Limit		*Plasma Protein Multicalibrator	
Blank response			
Difference Limit			
SD			

TEST PARAMETERS

No.	
Test	C4
Method	Endpoint
Direction	Ascend
Unit	mg/dL
Decimals	1
Prim. Wave.	340
Sec. Wave.	0
Sample Vol. (µL)	7.5
R1 Vol. (µL)	300
R2 Vol. (µL)	0
Line. Limit	
Incubation	0
Reaction	0 10
Antigen	<input type="checkbox"/> Check Antigen
Substrate	0
Response	L Limit 0 U Limit 0
S. Volume	0
Ratio	0
R1 Blank	L Limit 0 U Limit 0
Mix. R Blank	L Limit 0 U Limit 0
Linearity	L Limit 0.5 U Limit 100
Full Name	C4
Print No.	

For other instrument software version please contact with LINEAR