IgA at

Turbidimetric method

REF CT31700 2 x 40 mL	REF CT31702 6 x 40 mL
CONTENTS R1.Reagent 2 x 40 mL	CONTENTS R1.Reagent 6 x 40 mL
For <i>in vitro</i> diagnostic use only	

SUMMARY

IgA at is a quantitative turbidimetric assay^{1,2} for the measurement of IgA in human serum or plasma.

Anti-human IgA antibodies form insoluble complexes when mixed with samples containing IgA. The scattering light of the immunocomplexes depends of the IgA concentration in the patient sample, and can be quantified by comparison from a calibrator of known IgA concentration.

REAGENTS

R1 IgA at. Goat antibodies anti-human IgA, tris buffer 20 mmol/L, pH 8.2.

Precautions: The reagent contains sodium azide 0.95 g/L. Avoid any contact with skin or mucous.

PREPARATION

R1. Ready to use.

Calibration curve:

Prepare two fold dilutions of the Calibrator (1/1, 1/2, 1/4, 1/8) using NaCl 9 g/L as diluent. Use NaCl 9 g/L as Calibrator 0. Calculate the concentration of each IgA Calibrator dividing the Calibrator concentration by the corresponding dilution factor.

STORAGE AND STABILITY

1.  The reagents will remain stable until the expiration date printed on the label. For optimal stability store **tightly closed** at 2-8°C. Do not use the reagents after the expiration date.
2. Presence of particles, turbidity and/or the absorbance of blank reagent > 0.3 at 340 nm are sign of deterioration.

SAMPLE COLLECTION

Fresh serum and EDTA or heparinized plasma.

IgA in serum or plasma is stable 7 days at 2-8°C or 3 months at -20°C.

Samples with presence of fibrin should be centrifuged before testing. Highly hemolyzed or lipemic samples are not suitable for testing.

INTERFERENCES

- Bilirubin (10 mg/dL) does not interfere. Hemoglobin (4 g/L), rheumatoid factors (200 IU/mL) and lipemia (1.25 g/L) may affect the results. Other substances may interfere⁵.

ADDITIONAL EQUIPMENT

- Analyzer LIDA.

AUTHOMATIC TECHNIQUE

For automatic assays the application method is included.

Any application to an instrument should be validated to demonstrate that results meet the performance characteristics of the method.

It is recommended to validate periodically the instrument.

CALIBRATION

Recalibrate weekly, when a new lot of reagent is used, when control recovery falls out of the expected range or when adjustments are made to the instrument. A reagent blank should be run daily before sample analysis.

REFERENCE VALUES

Adults³: 70 – 400 mg/dL

Newborn⁴: 0 – 2.2 mg/dL

It is recommended that each laboratory establishes its own reference range.

QUALITY CONTROL

To ensure adequate quality control (QC), each run should include a set of controls (normal and abnormal) with assayed values handled as unknown. Each laboratory should establish its own Quality Control scheme and corrective actions if controls do not meet the acceptable tolerances.

DIAGNOSTIC CHARACTERISTICS³⁻⁵

Quantification of immunoglobulins in serum is important for the diagnosis of primary or secondary immunodeficiency, for monitoring of immunoglobulin therapy, and for following the clinical course of multiple myeloma⁴.

IgA represents about 10 to 15 % of the serum immunoglobulins. Its monomeric form has a similar structure that IgG, although between 10 and 15% of the total IgA is dimeric, more resistant to destruction by some pathogenic bacteria⁴.

The most important form of IgA is the secretory IgA, found in tears, sweat, saliva, milk, colostrums and gastrointestinal and bronchial secretions.

Congenital and acquired immunodeficiency are the most important causes of IgA deficit.

Polyclonal hyperimmunoglobulinemia (normal response to infections), skin, gut, respiratory and renal infections as well as cirrhosis are causes that increase IgA concentration⁴.

Increments of monoclonal IgA (paraprotein) are found in multiple myeloma, and proliferative disorders of plasma cells⁴.

PERFORMANCE CHARACTERISTICS

Performance characteristics are available on request.

NOTES

1. The linearity limit depends on the sample/reagent ratio, as well as the analyzer used. It will be higher by decreasing the sample volume, although the sensitivity of the test will be proportionally decreased.
2. Clinical diagnosis should not be made on findings of a single test result, but should integrate both clinical and laboratory data.

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REFERENCIAS

1. Narayanan S. *Clin Chem* 128: 1528-1531 (1982).
2. Price CP et al. *Ann Clin Biochem* 20: 1-14 (1983).
3. Dati F et al. *Eur J Clin Chem Clin Biochem* 34:517-520 (1996).
4. Tietz *Textbook of Clinical Chemistry*, 3rd Ed. Burtis CA, Ashwood ER. WB Saunders Co., (1999).
5. Young DS. *Effects of drugs on clinical laboratory tests*. 3th ed. AACC Press (1997).
6. Friedman and Young. *Effects of the disease on clinical laboratory tests*, 3th ed. AACC Press, 1997.

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CALIBRATION PARAMETERS

Rule:	Logistic 5 P	Calibrator 1:	Deionized Water
Factor K		Calibrator 2:	*PPM1
Replicates	2	Calibrator 3:	*PPM2
Interval		Calibrator 4:	*PPM3
Sensitivity		Calibrator 5:	*PPM4
Error Limit		*Plasma Protein Multicalibrator	
Blank response			
Difference Limit			
SD			

TEST PARAMETERS

No.	
Test	IgA
Method	Endpoint
Direction	Ascend
Unit	mg/dL
Decimals	1
Prim. Wave.	340
Sec. Wave.	0
Sample Vol. (µL)	3
R1 Vol. (µL)	450
R2 Vol. (µL)	0
Line. Limit	
Incubation	
Reaction	0 10
Antigen	<input type="checkbox"/> Check Antigen
Substrate	0
Response	L Limit 0 U Limit 0
S. Volume	
Ratio	
R1 Blank	L Limit 0 U Limit 0
Mix. R Blank	L Limit U Limit
Linearity	L Limit 4 U Limit 800
Full Name	IgA
Print No.	

For other instrument software version please contact with LINEAR