usCRP - Turbidimetric

Latex Turbidimetry

REF CT32500 1 x 50 mL CONTENTS R1.Reagent 1 x 40 mL R2.Reagent 1 x 10 mL
For <i>in vitro</i> diagnostic use only

PRINCIPLE

The latex particles coated with anti-CRP are agglutinated when they react with samples that contain low levels of C-reactive protein (CRP). The latex particles agglutination is proportional to the concentration of the CRP in the sample and can be measured by turbidimetry.¹

REAGENTS

R1. Diluent. Tris buffer, 20 mmol/L, pH 8.2.

R2. Latex. Latex particles coated with goat anti-human CRP, pH 8.0.

Precautions: The reagents contain sodium azide 0.95 g/L. Avoid any contact with skin or mucous.

PREPARATION

R1. Ready to use.

R2. Ready to use. Shake gently the vial before use.

Calibration curve: Prepare two fold dilutions of the Calibrator (1/1, 1/2, 1/4, 1/8) using NaCl 9 g/L as diluent. Use NaCl 9 g/L as Calibrator 0. Calculate the concentration of each usCRP Calibrator dividing the Calibrator concentration by the corresponding dilution factor.

STORAGE AND STABILITY

-  The reagents will remain stable until the expiration date printed on the label. For optimal stability stored **tightly closed** at 2-8°C. Do not use the reagents after the expiration date.
- Reagent deterioration: Presence of particles and turbidity.

SAMPLE COLLECTION

Fresh serum. Stable 7 days at 2-8°C or 3 months at -20°C. Samples with presence of fibrin should be centrifuged before testing. Hemolyzed or contaminated samples are not suitable for testing.

INTERFERENCES

- Bilirubin (40 mg/dL), hemoglobin (6 g/L), and rheumatoid factors (150 IU/mL), do not interfere. Lipemia (10 g/L, intralipid), interferes. Other substances may interfere⁸.

ADDITIONAL EQUIPMENT

- LIDA Analyzer.
- Laboratory equipment.
- iso-Clean Solution. Ref. CT18002.
- Linear us-CRP Calibrator Ref. 3925005.

AUTHOMATIC TECHNIQUE

For automatic assays the application is included. Any application to an instrument should be validated to demonstrate that results meet the performance characteristics of the method. It is recommended to validate periodically the instrument.

CALIBRATION

Recalibrate weekly, when a new lot of reagent is used, when control recovery falls out of the expected range or when adjustments are made to the instrument. A reagent blank should be run daily before sample analysis.

EXPECTED VALUES^{3,7,8}

Normal range in adults: < 5 mg/L.
 Cardiovascular risk in adults: < 3 mg/L.
 Each laboratory should establish its own reference range.

QUALITY CONTROL

To ensure adequate quality control (QC), each run should include a set of controls (normal and abnormal) with assayed values handled as unknowns. Each laboratory should establish its own Quality Control scheme and corrective actions if controls do not meet the acceptable tolerances.

CLINICAL SIGNIFICANCE²⁻⁴

La CRP is an acute-phase protein present in normal serum, which increases significantly after most forms of tissue injuries, bacterial and virus infections, inflammation and malignant neoplasia. During tissue necrosis and inflammation resulting from microbial infections, the CRP concentration can rise up to 300 mg/L in 12-24 hours. Elevations of CRP levels within the conventional reference range have been reported in several studies, where demonstrate that CRP may be also useful in detecting atherosclerotic process and providing important prognostic information about patients with asymptomatic heart disease, unstable angina, and myocardial infarction. More studies in apparently healthy people show that CRP concentration in serum rise long before traditional symptoms of heart and vascular diseases are noticed.

PERFORMANCE CHARACTERISTICS

Performance characteristics are available on request.

NOTES

- The linearity limit depends on the sample/reagent ratio, as well as the analyzer used. It will be higher by decreasing the sample volume, although the sensitivity of the test will be proportionally decreased.
- The patient should be metabolically stable before to proceed the assay. Any suspicion that the CRP may be greater than 30 mg/L, dilute sample 1/10 in CLNa 9 g/L before to run the assay.
- Clinical diagnosis should not be made on findings of a single test result, but should integrate both clinical and laboratory data.

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