HbA1c - Turbidimetric

Latex Turbidimetry

REF CT33000 1 x 45 mL CONTENTS R1.Reagent 1 x 30 mL R2a. Reagent 1 x 10 mL R2b. Reagent 1 x 5 mL	REF 3156005 1 x 125 mL CONTENTS R1.Reagent 1 x 125 mL
For <i>in vitro</i> diagnostic use only	

SUMMARY

This method directly determinates the hemoglobin A1c (HbA1c) in whole blood, using an antigen and antibody reaction. Total hemoglobin and HbA1c compete for the unspecific absorption rate to the latex particles (R1). When anti-human HbA1c monoclonal antibody is added (R2), latex-HbA1c-anti-human HbA1c antibody complex is formed. The presence of goat anti-mouse IgG polyclonal antibody causes the agglutination of the particles (complexes). The amount of agglutination is proportional to the concentration of the HbA1c in the sample and can be measured by turbidimetry.

REAGENTS

R1. Latex Suspension. Polystyrene particles in a buffer, stabilizer. pH 8.0 +/- 0.2

R2a. Antibody A. Mouse anti-human HbA1c monoclonal antibody in a buffer.

R2b. Antibody B. Goat anti-mouse IgG polyclonal antibody in a buffer.

Lyse Reagent. Ref: 3156005 Buffer.

Precautions: The reagents contain sodium azide 0.95 g/L. Avoid any contact with skin or mucous.

PREPARATION

R1. Ready to use. Mix thoroughly before use.

R2. Pour the total volume of R2b vial into the R2a vial. Stable 30 days at 2-8°C.

HbA1c Calibrator Ref: 3955005. **Calibration curve:** Use the 4 levels of HbA1c Calibrator. Use distilled water as calibrator 0.

STORAGE AND STABILITY

- The reagents will remain stable until the expiration date printed on the label. Opened vials are stable 1 month when stored at 2-8°C. For optimal stability store **tightly closed** at 2-8°C. Do not use the reagents after the expiration date. Do not interchange reagents from different lots.
- Reagent deterioration: Presence of particles and turbidity.

SAMPLE COLLECTION

Venous blood collected with EDTA (Note 1).

- Dispense 1 mL *Lyse Reagent* into labeled tubes (Note 2): Controls, or patient samples.
- Place 20 µL of well mixed blood (patient sample or controls) to each tube. Mix.
- Incubate for 5 minutes at room temperature until complete lysis is evident. Stable 10 days at 2-8°C.

INTERFERENCES

- Bilirubin (50 mg/dL), ascorbic ac. (50 mg/dL), triglycerides (2000 mg/dL), carbamylated Hb (7.5 mmol/L), acetylated Hb (5.0 mmol/L), uremia and labile intermediates (Schiff base), do not interfere. Hemoglobin variants HbA2, HbC and HbS do not interfere. Other rare hemoglobins (e.g. HbE) have not been studied
- Other substances or drugs may interfere.¹¹

ADDITIONAL EQUIPMENT

- LIDA Analyzer.
- Laboratory equipment.
- iso-Clean Solution. Ref. CT18002.
- Linear HbA1c Calibrator Ref. 3955005.

AUTHOMATIC TECHNIQUE

For automatic assays the application method is included. Any application to an instrument should be validated to demonstrate that results meet the performance characteristics of the method. It is recommended to validate periodically the instrument.

CALIBRATION

Recalibrate weekly, when a new lot of reagent is used, when control recovery falls out of the expected range or when adjustments are made to the instrument. A reagent blank should be run daily before sample analysis.

EXPECTED VALUES

Table below indicates the cut-off values established by the *Diabetes Control and Complications Trial Research Group* (DCCT) and have been adopted for a reference population (non-diabetic) and for the evaluation of the degree blood glucose control in diabetic patients^{2,5}.

DCCT / NGSP	IFCC	Degree of control
4.0 – 6.0	2.0 – 4.2	Non Diabetic
6.0 – 6.5	4.2 – 4.8	Goal
6.5 – 8.0	4.8 – 6.4	Good Control
> 8.0	> 6.4	Action suggested

Each laboratory should establish its own expected values (Note 3 and 4).

QUALITY CONTROL

To ensure adequate quality control (QC), each run should include a set of HbA1c controls (ref: 3955010.) with assayed values handled as unknown. Each laboratory should establish its own Quality Control scheme and corrective actions if controls do not meet the acceptable tolerances.

DIAGNOSTIC CHARACTERISTICS

HbA1c is the product of the irreversible condensation of glucose with the N-terminal residue of β-chain of hemoglobin A. The concentration of HbA1c is directly proportional to the concentration of glucose in the blood, and represents the integrated values for glucose over the preceding 6 to 8 week. HbA1c values are free of day-to day glucose fluctuations and unaffected by exercise or recent food ingestion. The determination of HbA1c is most commonly performed for the evaluation of glycemic control in diabetic patients, although is not adequate for diabetes mellitus control. A higher HbA1c value indicates poorer glycemic control. This process, reflects the average exposure of hemoglobin to glucose over an extended period. HbA1c in diabetic subjects use to be elevated 2-3 fold over the levels found in normal individuals¹. HbA1c serve as an indicator of metabolic control of the diabetic, since HbA1c levels approach normal values for diabetics in metabolic control.²⁻⁴

PERFORMANCE CHARACTERISTICS

Performance characteristics are available on request.

NOTES

- HbA1c in whole blood collected with EDTA is stable for one week at 2-8°C⁵.
- Plastic or glass tubes are acceptable.
- In using HbA1c to monitor diabetic patients, results should be interpreted individually. That is, the patient should be monitored against himself. There is a 3-4 week time lag before HbA1c reflects changes in blood glucose level.
- Clinical diagnosis should not be made on findings of a single test result, but should integrate both clinical and laboratory data.
- The linearity limit depends on the sample/reagent ratio, as well as the analyzer used. It will be higher by decreasing the sample volume, although the sensitivity of the test will be proportionally decreased.

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Linear HbA1c Calibrator. Code 3955005

CONTROL CALIDAD
QUALITY CONTROL

Linear HbA1c Control Level-I and Level II. Code 3955010

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CALIBRATION PARAMETERS

Rule:	Spline	Calibrator 1:	Deionized Water
Factor K		Calibrator 2:	HbA1c CAL 1 (*)
Replicates	2	Calibrator 3:	HbA1c CAL 2 (*)
Interval		Calibrator 4:	HbA1c CAL 3 (*)
Sensitivity		Calibrator 5:	HbA1c CAL 4 (*)
Coefficient		(*) Insert concentration	
Error Limit			
Blank response			
Difference Limit			
SD			

TEST PARAMETERS

No.		
Test	HbA1c	
Method	Endpoint	
Direction	Ascend	
Unit	%	
Decimals	1	
Prim. Wave.	670	
Sec. Wave.	0	
Sample Vol. (µL)	6	
R1 Vol. (µL)	200	
R2 Vol. (µL)	100	
Line. Limit	0	
Incubation	25	
Reaction	0	25
Antigen	<input type="checkbox"/> Check Antigen	
Substrate		
Response	L Limit 0	U Limit 0
S. Volume	0	
Ratio	0	
R1 Blank	L Limit 0	U Limit 0
Mix. R Blank	L Limit 0	U Limit 0
Linearity	L Limit 2	U Limit 16
Full Name	HbA1c	
Print No.		

For other instrument software version please contact with LINEAR