CREATININE E

Colorimetric enzymatic method


END POINT

REF KR10150 5 x 45 mL CONTENTS R1. Reagent 5 x 30 mL R2. Reagent 5 x 15 mL	REF KR10152 11 x 45 mL CONTENTS R1. Reagent 11 x 30 mL R2. Reagent 11 x 15 mL
For <i>in vitro</i> diagnostic use only	


SUMMARY

This procedure¹ is based upon an improved enzymatic method originally devised to evaluate creatinine in serum and urine. The test is performed in two steps.

In the first step creatine is removed during the first minutes of the preincubation stage of the sample with creatinase. In the second step the addition of creatininase, acting also as the reaction starter, hydrolyzes the creatinine of the sample in the presence of sarcosine oxidase (Sar OD) with production of hydrogen peroxide.

The hydrogen peroxide produced from the oxidase reaction is quantified by a Trinder's type reaction in which the chromogenic derivative HTIB* and 4-aminoantipyrine (4-AA) are condensed in the presence of peroxidase (POD) to form a red quinoneimine dye.

The rate of color development is proportional to the concentration of creatinine in the sample.

*Hydroxy-triiodobenzoic acid.

REAGENTS

R1. Chromogen reagent. TAPS buffer 35 mmol/L pH 8.0, creatinase \geq 300 μ kat/L, sarcosine oxidase \geq 100 μ kat/L, ascorbate oxidase \geq 20 μ kat/L, HTIB 2 mmol/L, detergent 0.5 g/L **R:36/37/38.**

R2. Enzyme reagent. TAPS buffer 50 mmol/L pH 8.0, creatininase \geq 400 μ kat/L, peroxidase \geq 15 μ kat/L, 4-amino antipyrine 2 mmol/L, potassium ferrocyanide 0.1 mol/L, detergent 1.5 g/L. **X_i R:36/37/38.**

PREPARATION

The Reagents are ready-to-use.

STORAGE AND STABILITY

Store at 2-8°C. The reagents are stable until the expiry date stated on the label. After daily use stored tightly closed and protected from light. On board the reagents are stable 15 days.

SAMPLE COLLECTION

Serum or heparinized plasma, and urine (see Notes).

In serum or plasma is stable up to 24 h at 2-8°C. Freeze for longer storage. In urine is stable up to 4 days at 2-8°C. Freeze for longer storage.

INTERFERENCES

- Bilirubin (>25 mg/dL), hemoglobin (0.5 g/dL), triglycerides (>2.0 g/dL), creatine (>20 mg/dL) and ascorbate (>30 mg/dL) may affect the results.
- A number of drugs are known to affect creatinine levels.²

INSTRUMENTATION AND MATERIALS

- KROMA analyzer.
- Laboratory equipment.
- Cleaning solutions Ref. KR18100, KR18200, KR18300, KR18400.
- Multicalibrator CC/H 10x3 mL Ref. CT19750.

AUTOMATED PROCEDURE

A graphic display pictures the specific sets corresponding to the technical application outlined for this test. Any new application should be validated to confirm that results meet the analytical performance of the method. It is recommended to validate periodically the instrument.

CALIBRATION

Recalibrate weekly, when a new lot of reagent is used, when control recovery falls out of the expected range or when adjustments are made to the instrument.

Two point calibration is recommended (S1: NaCl 9 g/L and S2: Calibrator). A reagent blank should be run daily before analysis.

RESULTS

Samples with concentrations higher than 20 mg/dL should be diluted 1:4 with saline and assayed again. Multiply the results by 4.

If results are to be expressed as SI units apply: mg/dL x 88.4 = μ mol/L

EXPECTED VALUES⁶

Serum, plasma

Men	0.70-1.20 mg/dL (62-106 μ mol/L)
Women	0.50-0.90 mg/dL (44-80 μ mol/L)

Urine

Men	14-26 mg/Kg/24-h (124-230 μ mol/Kg/24-h)
Women	11-20 mg/Kg/24-h (97-117 μ mol/Kg/24-h)

It is recommended that each laboratory establishes its own reference range.

QUALITY CONTROL

To ensure adequate quality control (QC), each run should include a set of controls (normal and abnormal) Ref. CT19800-CT19850 with assayed values handled as unknowns.

Each laboratory should establish its own Quality Control scheme and corrective actions if controls do not meet the acceptable tolerances.

DIAGNOSTIC CHARACTERISTICS

Elevated levels of creatinine in serum are usually associated with renal diseases, especially those related to GFR such as glomerular nephritis. Therefore, the clinical significance of the creatinine level in plasma or serum is usually determined in conjugation with the plasma urea level since there is an increase in both levels in postrenal azotemia, while the CC, or urine levels, are diminished.

Clinical diagnosis should not be made on findings of a single test result, but should integrate both clinical and laboratory data.

NOTES

1. Creatinine in urine may be assayed on fresh random samples. No special preparation of the patient is needed. Dilute specimen 1:50 with distilled water before the assay. Multiply the result by 50.

PERFORMANCE CHARACTERISTICS

Performance characteristics are available on request.

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