# CREATINE KINASE BR

CK NAC-ACTIVATED


UV enzymatic method

KINETIC

REF KR10160	REF KR10162
1 x 40 mL	5 x 40 mL
<b>CONTENTS</b>	<b>CONTENTS</b>
R1. Reagent 1 x 32 mL R2. Reagent 1 x 8 mL	R1. Reagent 5 x 32 mL R2. Reagent 5 x 8 mL
For <i>in vitro</i> diagnostic use only	

## SUMMARY

Creatine kinase (CK) catalyzes the reaction between creatine phosphate (CP) and adenosine 5'-diphosphate (ADP) with formation of creatine and adenosine 5'-triphosphate (ATP). The latter phosphorylates glucose to glucose-6-phosphate (G6P) in the presence of hexokinase (HK). G6P is oxidized to Gluconate-6P in the presence of reduced nicotinamide-adenine dinucleotide phosphate (NADP) in a reaction catalyzed by glucose-6-phosphate dehydrogenase (G6P-DH). The conversion is monitored kinetically at 340 nm by the rate of increase in absorbance resulting from the reduction of NADP to NADPH proportional to the activity of CK present in the sample. In this test<sup>1,2</sup> the presence of N-acetylcysteine (NAC) allows the optimal activation of the enzyme.

This test has been formulated according to the standardized method described by IFCC. Clin Chem Lab Med 2002; 40(6) : 635-642.

## REAGENTS

**R1. Buffer/Glucose/NAC.** Imidazol buffer 100 mmol/L pH 6.7, glucose 20 mmol/L, NAC 20 mmol/L, magnesium acetate 10 mmol/L, NADP 2.5 mmol/L, HK  $\geq 4$  KU/L, EDTA 2 mmol/L.

**R2. Substrate/Coenzymes.** CP 30 mmol/L, AMP 5 mmol/L, ADP 2 mmol/L, di(adenosine-5') pentaphosphate 10  $\mu$ mol/L, G6P-DH  $\geq 1.5$  KU/L.

## PREPARATION

The Reagents are ready-to-use.

## STORAGE AND STABILITY

 Store at 2-8°C. The Reagents are stable until the expiry date stated on the label. After daily use stored tightly closed and protected from light. On the board the reagents are stable 19 days. Discard the reagent if the blank presents an absorbance over 0.400 at 340 nm against distilled water.

## SAMPLE COLLECTION

Serum. Stable 8 days at 2-8°C or 1 month at -20°C. Chill the samples as rapidly as possible after collection. Moderately or severely hemolyzed specimens are unsatisfactory for testing, as well as plasmas containing EDTA, heparin, citrate, or fluoride since they may produce unpredictable reaction rates.

## INTERFERENCES

- Lipemia (intralipid >5 g/L) may affect the results.
- Bilirubin (<20 mg/dL), hemoglobin (<10 g/L), do not interfere.
- Other drugs and substances may interfere<sup>3</sup>.

## INSTRUMENTATION AND MATERIALS

- KROMA analyzer.
- Laboratory equipment.
- Cleaning solutions Ref. KR18100, KR18200, KR18300, KR18400.

## AUTOMATED PROCEDURE

A graphic display pictures the specific sets corresponding to the technical application outlined for this test. Any new application, to the instrument should be validated to confirm that results meet the analytical performance of the method.

It is recommended to validate periodically the instrument.

## CALIBRATION

No calibration is required. Activity is calculated against theoretical factor.

## CALCULATION

Patient values are calculated automatically by theoretical factor.

## RESULTS

Samples with  $\Delta A/\text{min}$  exceeding 0.270 a 340 nm should be diluted 1:10 with saline and assayed again. Multiply the results by 10.

If results are to be expressed as SI units apply: U/L x 16.67 = nkat/L

## EXPECTED VALUES<sup>1</sup>

Serum

Temperature	25°C	30°C	37°C
Men	$\leq 65$ U/L (1083 nkat/L)	$\leq 105$ U/L (1750 nkat/L)	$\leq 174$ U/L (2900 nkat/L)
Women	$\leq 55$ U/L (917 nkat/L)	$\leq 80$ U/L (1334 nkat/L)	$\leq 140$ U/L (2334 nkat/L)
Children	$\leq 94$ U/L (1570 nkat/L)	$\leq 150$ U/L (2500 nkat/L)	$\leq 225$ U/L (3750 nkat/L)

It is recommended that each laboratory establishes its own reference range.

## QUALITY CONTROL

To ensure adequate quality control (QC), each run should include a set of controls (normal and abnormal) Ref. CT19800-CT19850 with assayed values handled as unknowns.

Each laboratory should establish its own Quality Control scheme and corrective actions if controls do not meet the acceptable tolerances.

## DIAGNOSTIC CHARACTERISTICS

Creatine kinase (CK) values are high in patients with myocardial infarction, progressive muscular dystrophy, alcoholic myopathy, and delirium tremens, but normal in patients with hepatitis and other forms of liver disease. The high values in patients with hypothyroidism reflect the muscle changes in this condition. Although CK is found almost exclusively in myocardium, muscle, and brain and early reports suggested it to be an almost specific index of injury of myocardium and muscle, more recent reports indicate that, inexplicably high serum CK values can occur in patients with pulmonary infarction and pulmonary edema. Specificity of CK assay is enhanced measurement of its isoenzymes.<sup>4-6</sup> Clinical diagnosis should not be made on findings of a single test result, but should integrate both clinical and laboratory data.

## PERFORMANCE CHARACTERISTICS

Performance characteristics are available on request.

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