# LDL-CHOLESTEROL

DIRECT

*Enzymatic colorimetric method*

ENDPOINT

<b>REF KR10290</b> 5 x 40 mL <b>CONTENTS</b> R1. Reagent 5 x 30 mL R2. Reagent 5 x 10 mL	<b>REF KR10292</b> 1 x 40 mL <b>CONTENTS</b> R1. Reagent 1 x 30 mL R2. Reagent 1 x 10 mL
For <i>in vitro</i> diagnostic use only	

## SUMMARY

The assay is a homogeneous method for directly measuring LDL-C concentrations in serum or plasma, without the need for any off-line pretreatment or centrifugation steps.

The method is in a two-reagent format and depends on the properties of a unique detergent. This detergent (Reagent 1) solubilizes only the non LDL lipoprotein particles. The cholesterol released is consumed by cholesterol esterase and cholesterol oxidase in a non color forming reaction. A second detergent (Reagent 2) solubilizes the remaining LDL particles and a chromogenic coupler allows for color formation. The enzyme reaction with LDL-C in the presence of the coupler produces color that is proportional to the amount of LDL cholesterol present in the sample. The intensity of the color formed is proportional to the LDLc concentration in the sample.

## REAGENTS

**R1. Reagent 1.** Buffer. Detergent 1 <1%, Cholesterol esterase (*Pseudomonas* sp.) <1500 U/L, Cholesterol oxidase (*Cellulomonas* sp.) <1500 U/L, Peroxidase (Horseradish) <1300 ppg U/L, 4-Aminoantipyrine <0.1%, Ascorbic oxidase (*Curcubita* sp.) <3000 U/L. Preservative.

**R2. Reagent 2.** Buffer (pH 6.3). Detergent 2 <1%, N,N-bis (4-sulfobutyl) - m-toluidine, disodium (DSBmT) <1.0 mM. Preservative.

**PREPARATION** The Reagents are ready to use.

## STORAGE AND STABILITY

Store at 2-8°C. Stable until the expiry date stated on the label. After daily use stored tightly closed and protected from light. Do not use reagents over the expiration date. Do not freeze. On board are stables 4 weeks. Discard if appear signs of deterioration:

- Presence of particles and turbidity.
- Inability to recover control values.

## SAMPLE COLLECTION

Serum, EDTA or heparinized plasma obtained by the patient after an overnight fast. Remove from cells within 3 hours of venipuncture. Samples may be kept at 2-8°C for 2 weeks or at -20°C for 3 months. Anticoagulants containing citrate should not be used.

## INTERFERENCES

- Ascorbic Acid 50 mg/dL
- Hemoglobin 500 mg/dL
- Bilirubin 20 mg/dL
- Gamma-Globulins 5000 mg/dL
- Other drugs and substances may interfere.

Concentration with no significant ( $\pm 10\%$ ) interference.

## INSTRUMENTATION AND MATERIALS

- KROMA analyzer.
- Laboratory equipment.
- Cleaning solutions Ref. KR18100, KR18200, KR18300, KR18400.
- HDL/LDL-Cholesterol Calibrator Ref. 1972005.

## AUTOMATED PROCEDURE

A graphic display pictures the specific sets corresponding to the technical application outlined for this test. Any new application, to the instrument should be validated to confirm that results meet the analytical performance of the method. It is recommended to validate periodically the instrument.

## CALIBRATION

Recalibrate weekly, when a new lot of reagent is used, when control recovery falls out of the expected range or when adjustments are made to the instrument. Two point calibration is recommended (S1: NaCl 9 g/L and S2: Calibrator). A reagent blank should be run daily before sample analysis.

## RESULTS

If results are to be expressed as SI units apply:  $\text{mg/dL} \times 0.0259 = \text{mmol/L}$

## EXPECTED VALUES

Risk group classification according with LDL-C levels.

LDL-Cholesterol levels	RISK
< 100 mg/dL (2.59 mmol/L)	Optimal
100-129 mg/dL (2.59-3.34 mmol/L)	Near optimal
130-159 mg/dL (3.37-4.12 mmol/L).	Borderline
160-189 mg/dL (4.14-4.89 mmol/L)	High
$\geq 190$ mg/dL (4.92 mmol/L)	Very high

It is recommended that each laboratory establishes its own reference range.

## QUALITY CONTROL

To ensure adequate quality control (QC), each run should include a set of controls (normal and abnormal) with assayed values handled as unknowns.

Each laboratory should establish its own Quality Control scheme and corrective actions if controls do not meet the acceptable tolerances.

## DIAGNOSTIC CHARACTERISTICS

Laboratory investigations, epidemiology and genetic factors of hypercholesterolemia indicate that elevated LDL cholesterol is a major cause of coronary heart disease (CHD). The relationship between LDL-cholesterol levels and CHD risk is continuous over a broad range of LDL levels from low to high.

Recent clinical trials robustly show that LDL-lowering therapy reduces risk for CHD. For these reasons the Adult Treatment Panel III continues to identify elevated LDL cholesterol as the primary target of cholesterol-lowering therapy.

Clinical diagnosis should not be made on findings of a single test result, but should integrate both clinical and laboratory data.

## PERFORMANCE CHARACTERISTICS

Performance characteristics are available on request.

## BIBLIOGRAPHY

1. Gotto AM, Lipoprotein metabolism and the etiology of hyperlipidemia, Hospital Practice, 23: Suppl.1, 4 (1988).
2. Crouse JR et al., Studies of low density lipoprotein molecular weight in human beings with coronary artery disease, J. Lipid Res.,26; 566 (1985).
3. Badimon JJ, Badimon L., Fuester V., Regression of Atherosclerotic Lesions by High-Density Lipoprotein Plasma Fraction in the Cholesterol-Fed Rabbit, Journal of Clinical Investigation, 85:1234-41 (1990).
4. Castelli WP et al., Cholesterol and other lipids in coronary heart disease, Circulation, 55; 767 (1977).
5. Barr DP, Russ EM, Eder HA, Protein-lipid relationships in human plasma, Am. J. Med.,11;480 (1951).

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## REFERENCIAS

- Gotto AM, Lipoprotein metabolism and the etiology of hyperlipidemia, Hospital Practice, 23: Suppl.1, 4 (1988).
- Crouse JR et al., Studies of low density lipoprotein molecular weight in human beings with coronary artery disease, J. Lipid Res.,26; 566 (1985).
- Badimon JJ, Badimon L., Fuester V., Regression of Atherosclerotic Lesions by High-Density Lipoprotein Plasma Fraction in the Cholesterol-Fed Rabbit, Journal of Clinical Investigation, 85:1234-41 (1990).
- Castelli WP et al., Cholesterol and other lipids in coronary heart disease, Circulation, 55; 767 (1977).
- Barr DP, Russ EM, Eder HA, Protein-lipid relationships in human plasma, Am. J. Med.,11;480 (1951).

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## KROMA

**General Info**

Name:  Code:  Barcode:  Unit:  Decimal digits:

Type:   No. of rgt:

Visible Mode:  Notes:

---

**Filters**

F1:  F2:

**Volumes [microlitres]**

Sample:

R1:  R2:  R3:

**Bottle sizes**

R1:  R2:  R3:

**Incubation/Reading time [sec]**

Substrate/Sample Start:

R1,S -> R2:  2nd Incub. time:  Final incub.:

**Reference range**

Sample type:

Patient type	Min	Max
Female	0.00	0.00
Male	0.00	0.00

**Instrument Factor (Y = aX + b)**

a:  b:

**Controls**

C1  C2  C3

No. of SD for QC ref range:

---

**Reagents**

Blank Abs (min; max):    Include blank in calc.

Reagent linearity:

Detection limit:

Dilution settings

Name	N. of Standards	N. of rep.	Units	Stability
LDL-Cholesterol Direct	2	3	mg/dL	

Dil.ratio	Std value	O.D.	Reagent ...
	*		

Factor	Factor Min	Factor Max

\* entered by the user

KR10290-1/0907

## KROMA PLUS

**General Info**

Name:  Code:  Barcode:  Unit:  Decimal digits:

Type:    Bichromatic No. of rgt:  Multiply pre-diluted result:

Visible Mode:  Notes:

**Filters**

F1:  F2:

**Incubation/Reading time [sec]**

Final incub.:

**Reference range**

Sample type:

Patient type	Min	Max
Female	0.00	0.00
Male	0.00	0.00
Paediatric	0.00	0.00

**Volumes [microlitres]**

Sample:

R1:  R1 Bottle sizes:

R2:  R2 Bottle sizes:

**Instrument Factor (Y = aX + b)**

a:  b:

**Controls**

C1  C2  C3

No. of SD for QC ref range:

**Reagents**

Blank Abs (min; max):    Include blank in calc.

Reagent linearity:

Detection limit:

Show reagent blank history

Dilution settings

Name	N. of Standards	N. of rep.	Units	Stability
LDL-Cholesterol Direct	2	3	mg/dL	

Dil.ratio	Std value	O.D.	Reagent ...
*			

Factor	Factor Min	Factor Max

\* entered by the user

KR10290-1/0907