

CERTIFICATE OF ANALYSIS	IgA id
--------------------------------	---------------

REF	3300305	1 x 15 test			

LOT F041.15
 2-8°C

2017-02
IVD

FUNCTIONALITY

REAGENT	TEST	SPECIFICATIONS	RESULTS
IgA id	Slope	--	0.141700151
	Interp. Coeff. r ²	0.97-1.00	0.99305842
	Interception	10.0 – 11.0	10.5923529

CONTROL SERA	SPECIFICATIONS <i>mg/dl</i>	RESULTS <i>mg/dl</i>
1	276 – 374	322.2
2	221 -299	251.6
3	166 – 224	196.5
4	110 – 150	138.7
5	55 - 75	61.9

This lot was tested by test methods described in the manufacturer's package insert.
 These results are applicable to this lot and sequential sublots, for which a capital letter have been assigned after the lot number.

This Certificate of Analysis is a computer printout and has therefore not been signed by hand.

This lot meets our specifications and is released for sale	APPROVED	Signed - Date of Manufacture
	(APPROVED / NOT APPROVED)	Alicia Abellan Quality Control Department 11 January 2017



QUALITY SYSTEM CERTIFIED
ISO 9001 ISO 13485



LINEAR CHEMICALS, S.L.U. Joaquim Costa 18 2ª planta. 08390 Montgat (Barcelona) SPAIN
Telf. (+34) 934 694 990; E-mail: info@linear.es ; website: www.linear.es