

CERTIFICATE OF ANALYSIS	IgM id
--------------------------------	---------------

REF

LOT F054.17
 2-8°C

2019-02
IVD

FUNCTIONALITY

REAGENT	TEST	SPECIFICATIONS	RESULTS
IgA id	Slope	--	0.195420788
	Interp. Coeff. r ²	0.97-1.00	0.99096844
	Interception	10.0 – 11.0	10.0023529

CONTROL SERA	SPECIFICATIONS <i>mg/dl</i>	RESULTS <i>mg/dl</i>
1	96-130	109.3
2	77-104	92.6
3	58-78	61.9
4	38-52	43.4
5	19-26	22.7

This lot was tested by test methods described in the manufacturer's package insert. These results are applicable to this lot and sequential sublots, for which a capital letter have been assigned after the lot number.

This Certificate of Analysis is a computer printout and has therefore not been signed by hand.

This lot meets our specifications and is released for sale	APPROVED	Signed - Date of Manufacture,
	(APPROVED / NOT APPROVED)	Alícia Abellan Quality Control Department 20 June 2017

