


CERTIFICATE OF ANALYSIS	IgG id
-------------------------	--------

REF	3300505	15 Test	x		

**LOT** D271.16  
 2-8°C

 2018-09  
**IVD**

### FUNCTIONALITY

REAGENT	TEST	SPECIFICATIONS	RESULTS
IgG id	Slope	--	0.0284211
	Interp. Coeff. $r^2$	0.97 - 1	0.99566251
	Interception	10.0 – 11.0	10.6784706

CONTROL SERA	SPECIFICATIONS <i>mg/dl</i>	RESULTS <i>mg/dl</i>
1	1382-1871	1603
2	1106-1497	1299
3	829-1123	1021
4	553-749	650
5	276-375	305

This lot was assayed by the test method described in the manufacturer's package insert.  
 These results are applicable to this lot and subsequent sublots, for which a capital letter has been assigned after the lot number.

This Certificate of Analysis is a computer printout and has therefore not been signed by hand.

This lot meets our specifications and is released for sale	APPROVED	Signed - Date of Manufacture,
	(APPROVED / NOT APPROVED)	Alícia Abellan Quality Control Department 12 January 2017

