

<b>CERTIFICATE OF ANALYSIS</b>	<b>TRANFERRIN id</b>
--------------------------------	----------------------

<b>REF</b>	3300805	1 x 15 test			

<table border="1" style="border-collapse: collapse;"> <tr> <td style="padding: 2px;"><b>LOT</b></td> <td style="padding: 2px;">G036.16</td> </tr> <tr> <td style="text-align: center;"></td> <td style="padding: 2px;">2-8°C</td> </tr> </table>	<b>LOT</b>	G036.16		2-8°C	<table style="border-collapse: collapse;"> <tr> <td style="text-align: center;"></td> <td style="padding: 2px;">2018.02</td> </tr> <tr> <td style="text-align: center;"><table border="1" style="border-collapse: collapse;"><tr><td style="padding: 2px;"><b>IVD</b></td></tr></table></td> <td></td> </tr> </table>		2018.02	<table border="1" style="border-collapse: collapse;"><tr><td style="padding: 2px;"><b>IVD</b></td></tr></table>	<b>IVD</b>	
<b>LOT</b>	G036.16									
	2-8°C									
	2018.02									
<table border="1" style="border-collapse: collapse;"><tr><td style="padding: 2px;"><b>IVD</b></td></tr></table>	<b>IVD</b>									
<b>IVD</b>										

**FUNCTIONALITY**

REAGENT	TEST	SPECIFICATIONS	RESULTS
<b>TRANFERRIN id</b>	Slope	--	0.183129645
	Interp. Coeff. r <sup>2</sup>	0.97-1.00	0.99394764
	Interception	10.0 – 11.0	10.5476471

This lot was tested by test methods described in the manufacturer's package insert.  
 These results are applicable to this lot and sequential sublots, for which a capital letter have been assigned after the lot number.

This Certificate of Analysis is a computer printout and has therefore not been signed by hand.

This lot meets our specifications and is released for sale	APPROVED	Signed - Date of Manufacture Mercè Oriol 04 April 2016 Quality Control Department
	(APPROVED / NOT APPROVED)	

