



CERTIFICATE OF ANALYSIS	C3 id
--------------------------------	--------------

REF	3339005	15 Test	x
-----	---------	---------	---

 LOT AJ354.16 2-8°C	 2018-12 IVD
--	---

FUNCTIONALITY

REAGENT	TEST	SPECIFICATIONS	RESULTS
C3 id	Slope	--	0.332491795
	Interp. Coeff. r^2	0.97 - 1	0.99285812
	Interception	10.0 – 11.0	10.9342353

CONTROL SERA	SPECIFICATIONS <i>mg/dl</i>	RESULTS <i>mg/dl</i>
1	198-268	227.2
2	158-214	189.6
3	119-161	150.1
4	79-107	94.2
5	40-54	48.4

This lot was assayed by the test method described in the manufacturer's package insert.

These results are applicable to this lot and subsequent sublots, for which a capital letter has been assigned after the lot number.

This Certificate of Analysis is a computer printout and has therefore not been signed by hand.

This lot meets our specifications and is released for sale	APPROVED	Signed - Date of Manufacture, Alicia Abellan Quality Control Department 20 June 2017
	(APPROVED / NOT APPROVED)	

