

CERTIFICATE OF ANALYSIS	STEL 3-Diluent 20 L
--------------------------------	----------------------------

REF	3802005	20 L	X
------------	---------	------	---

Cont	DIL	X	LOT	38036		2020.04	IVD
-------------	------------	----------	------------	-------	-----------------------------------------------------------------------------------	---------	------------

PHYSICAL AND CHEMICAL CHARACTERISTICS

TEST		SPECIFICATIONS	RESULT	
Appearance	DIL	Transparent liquid; No precipitation, particles, or floccule	ACCEPTED	X
pH	DIL	7.5 ± 0.2	ACCEPTED	X

(ACCEPTED / N/A)

FUNCTIONALITY

INSTRUMENT	STEL 3
-------------------	---------------

TEST	SPECIFICATIONS	RESULTS
Blank counting	Result of a measurement using Hematology Analyzer, test value should be: White blood cell count(WBC) ≤ 0.3*10 ⁹ /L, Red blood cell count(RBC) ≤ 0.05*10 ¹² /L, Platelet count (PLT) ≤ 10*10 ⁹ /L, Hemoglobinometry (HGB) ≤ 2g/L.	WBC:0 0 0 RBC:0 0 0 PLT:0 0 0 HGB:0 0 0
Osmolality Value	Osmolality value should not exceed ±10 mOsm/kg of nominal value	307
Conductivity Value	Conductivity value should not exceed ±0.5 ms/cm of nominal value	16.94
Component	NaCl < 0.6%, Stabilizer < 1.0%	pass

(ACCEPTED / N/A)

This lot was assayed by the test method described in the manufacturer's package insert. These results are applicable to this lot and subsequent sublots, for which a capital letter has been assigned after the lot number.

This Certificate of Analysis is a computer printout and has therefore not been signed by hand.

This lot meets our specifications and is released for sale	APPROVED	Signed Date of Manufacture Mercè Oriol Quality Control Dept. 02 May 2017
-------------------------------------------------------------------	----------	----------------------------------------------------------------------------------------------------

