

<b>CERTIFICATE OF ANALYSIS</b>	<b>STEL 3- Extra Cleaner 100 mL</b>
--------------------------------	-------------------------------------

<b>REF</b>	3802020	100 mL	X
------------	---------	--------	---

<b>Cont</b>	Extra Cleaner	<b>X</b>	<b>LOT</b>	38037		2020/05	<b>IVD</b>
-------------	---------------	----------	------------	-------	---	---------	------------

**PHYSICAL AND CHEMICAL CHARACTERISTICS**

TEST	SPECIFICATIONS	RESULT	
Appearance	Transparent liquid; No precipitation, particles, or floccule	ACCEPTED	X
pH	7 ± 0.5	ACCEPTED	X

(ACCEPTED / N/A)

**FUNCTIONALITY**

<b>INSTRUMENT</b>	<b>STEL 3</b>
-------------------	---------------

TEST	SPECIFICATIONS	RESULTS
Blank counting	Result of a measurement using Hematology Analyzer, test value should be: White blood cell count(WBC)≤0.3*10 <sup>9</sup> /L, Red blood cell count(RBC) ≤0.05*10 <sup>12</sup> /L, Platelet count (PLT) ≤10*10 <sup>9</sup> /L, Hemoglobinometry (HGB) ≤2g/L.	WBC:0 0 0 RBC:0 0 0 PLT:0 0 0 HGB:0 0 0
Component	Buffer solution<0.3%,Protease<0.2%,	pass

(ACCEPTED / N/A)

This lot was assayed by the test method described in the manufacturer's package insert. These results are applicable to this lot and subsequent sublots, for which a capital letter has been assigned after the lot number.

This Certificate of Analysis is a computer printout and has therefore not been signed by hand.

<b>This lot meets our specifications and is released for sale</b>	APPROVED	<b>Signed</b> <b>Date of Manufacture</b> Mercè Oriol Quality Control Dept. 07 June 2017
---	----------	---

