



CERTIFICATE OF ANALYSIS	B-hCG strip
--------------------------------	--------------------

REF	4120050	50 Test	X
-----	---------	---------	---

Cont	Strips	LOT	F703025	 2019-02
			2-30 °C	IVD

PHYSICAL CHARACTERISTICS

REAGENT	TEST	SPECIFICATIONS	RESULTS
STRIP	APPEARANCE	COMPACT STRIP	PASS

(PASS / FAILURE)

FUNCTIONALITY TEST

ANALYSIS		STANDARDS	SPECIFICATIONS	RESULTS
Specificity:	300 mIU/ml	LH	Negative	Negative
	100 mIU/ml	hFSH	Negative	Negative
	1000 µIU/ml	hTSH	Negative	Negative
	0 mIU/ml	hCG (4 TH IPR)	Negative	Negative
Sensitivity claimed:	20 mIU/ml	hCG (4 TH IPR)	Positive	Positive
Prozone:	600 IU/ml	hCG (4 TH IPR)	Positive	Positive
Pass / No Pass	PASS			

This lot was tested by test methods described in the manufacturer's package insert.

These results are applicable to this lot and sequential sublots, for which a capital letter have been assigned after the lot number.

This Certificate of Analysis is a computer printout and has therefore not been signed by hand.

This lot meets our specifications and is released for sale	APPROVED	Signed Date of Manufacture, Mercè Oriol Quality Control Dept. 02 May 2017
	(APPROVED / NOT APPROVED)	

