

CERTIFICATE OF ANALYSIS	FECAL OCCULT BLOOD CASSETTE
--------------------------------	------------------------------------

REF

4325225	25 Test		
---------	---------	--	--

Cont FOB test device
Specimen collection tubes with extraction buffer

LOT T703026

 2019-02

 2-30 °C

IVD

PHYSICAL CHARACTERISTICS

REAGENT	TEST	SPECIFICATIONS	RESULTS
DEVICE	APPEARANCE	COMPACT DEVICE	PASS

(PASS / FAILURE)

FUNCTIONALITY TEST

ANALYSIS	SPECIFICATIONS	TEST RESULTS
SPECIFICITY	Extraction buffer	Negative
SENSITIVITY	40 ng/ml hHb	Positive
PROZONE	0.5 mg/ml hHb	Positive
RESULTS	PASS	

(PASS / FAILURE)

This lot was tested by test methods described in the manufacturer's package insert.

These results are applicable to this lot and sequential sublots, for which a capital letter have been assigned after the lot number.

This Certificate of Analysis is a computer printout and has therefore not been signed by hand.

This lot meets our specifications and is released for sale	APPROVED	Signed Date of Manufacture, Mercè Oriol Quality Control Dept. 29 March 2017
	(APPROVED / NOT APPROVED)	

