

CERTIFICATE OF QUALITY CONTROL

PLUS-SED auto (Tapval)	REF	5220060
-------------------------------	------------	----------------

LOT	180062		2019-05
------------	---------------	---	----------------



LOT INSPECTION RESULTS

TEST	SPECIFICATION	RESULT			
Citrate 3Na solution (flame photometer [Na ⁺])	0.1-0.136 Mol solution/l (0.3-0.4 Mol [Na ⁺]/l) Aqueous solution equivalent to 3.8% (w / v)	PASS	X	FAIL	
Sample volume	1.5 mL	PASS	X	FAIL	
Tube (gravimetric) ($\delta \approx 1.012$)	0.450-0.550 (0.455-0.555 gr) Citrate 3Na solution / tube	PASS	X	FAIL	
Visual inspection	UNE-66020-1 (\approx ISO 2859-1) NII-NCA/AQL:0. 1-0.4-1	PASS	X	FAIL	
Functionality	Level I (1 h): 7-21 Level 2 (1 h): 68-118	PASS	X	FAIL	

This lot was assayed under test method described in Linear Chemicals internal procedures

This Certificate of Analysis is a computer printout and has therefore not been signed by hand.

This lot meets our specifications and is released for sale	APPROVED	Signed Date of Manufacture, Mercè Oriol 12/12/2017 Quality Control
	(APPROVED / NOT APPROVED)	

