

## CERTIFICATE OF ANALYSIS

## Free T3 EIA 96 Test

## REF

6107211

96 Test

Cont	ft3 Ab Coated Wells	1
	Standard Liquid	1 set
	Conjugate	10,5 ml
	Reagent A	13 ml
	Reagent B	13 ml
	Stop Solution	10 ml

LOT

RN-57436



2018-12



2-8°C

IVD

## STANDARD CURVE

Concentration (ng/mL)	O.D 450 nm
0.0	2.490
1,8	2.182
3,3	1.777
5,7	1.190
8,6	0.761
18	0.396

This standard Curve was obtained on the date of performance testing of the above components. It is for the purpose of illustration only. Each laboratory should establish its own standard curve

## CONTROL VALUES

Control No.	Mean (ng/mL)	Range
1	2.45	1.76-3.50
2	4.83	3.31-6.57
3	7.25	4.72-9.38

This lot was tested by test methods described in the manufacturer's package insert.

These results are applicable to this lot and sequential sublots, for which a capital letter have been assigned after the lot number.

This Certificate of Analysis is a computer printout and has therefore not been signed by hand.

This lot meets our specifications and is released for sale	APPROVED	Signed – Date of manufacture Alícia Abellan Quality Control Dept. 30/08/2017
	(APPROVED / NOT APPROVED)	

