

CERTIFICATE OF ANALYSIS	T4 EIA 96 Test
--------------------------------	-----------------------

REF

6107215	96 Test	X
---------	---------	---

Cont	T4 Ab Coated Wells	1
	Standard Liquid	1
	Conjugate concentrate	1.3 ml
	Conjugate Diluent	13 ml
	TMB reagent	11 ml
	Stop Solution	11 ml

LOT RN-61141 2018-12
 2-8°C **IVD**

STANDARD CURVE

Concentration (ug/dL)	O.D 450 nm
0	2.669
2	1.687
5	1.188
10	0.832
15	0.576
25	0.307

This curve was obtained on the date of performance testing of the components. It is for the only purpose of illustration. Each laboratory should establish its own standard curve

CONTROL VALUES

Control No.	Mean (ug/dL)	Range
1	3.89	3.08-5.03
2	8.17	6.39-10.9
3	11.2	9.27-15.6

This lot was tested by test methods described in the manufacturer's package insert.

These results are applicable to this lot and sequential sublots, for which a capital letter have been assigned after the lot number.

This Certificate of Analysis is a computer printout and has therefore not been signed by hand.

This lot meets our specifications and is released for sale	APPROVED	Signed Date of Manufacture, Alícia Abellan 30/08/2017 Quality Control Dept.
	(APPROVED / NOT APPROVED)	

