

CERTIFICATE OF ANALYSIS	Free T4 EIA 96 Test
--------------------------------	----------------------------

REF

6107221	96 Test		
---------	---------	--	--

Cont	fT4 Ab Coated Wells	1
	Standard Liquid	1
	Enzyme Conjugate	10.5 ml
	Color Reagent A	13 ml
	Color Reagent B	13 ml
	Stop Solution	10 ml

LOT RN-57905 2019-04
 2-8°C **IVD**

STANDARD CURVE

Concentration (ng/dL)	O.D 450 nm
0.00	2.557
0.35	2.274
0.95	1.894
1.9	1.257
3.5	0.797
7.5	0.291

This curve was obtained on the date of performance testing of the components. It is for the only purpose of illustration. Each laboratory should establish its own standard curve

CONTROL VALUES

Control No.	Mean (ug/dL)	Range
1	0.555	0.139-1.04
2	2.18	1.23-2.65
3	4.92	2.42-5.40

This lot was tested by test methods described in the manufacturer's package insert.

These results are applicable to this lot and sequential sublots, for which a capital letter have been assigned after the lot number.

This Certificate of Analysis is a computer printout and has therefore not been signed by hand.

This lot meets our specifications and is released for sale	APPROVED	Signed Date of Manufacture, Alícia Abellan Quality Control Dept. 03 January 2018
	(APPROVED / NOT APPROVED)	

