

CERTIFICATE OF ANALYSIS	PROLACTIN EIA 96 Test
--------------------------------	------------------------------

REF

6107305	96 Test	X
---------	---------	---

Cont	Prolactin Ab Coated Wells	1
	Standard Lyophilized	1
	Conjugate reagent	13 ml
	TMB reagent	11 ml
	Stop Solution	11 ml

LOT RN-57157 2018-11
 2-8°C **IVD**

STANDARD CURVE

Concentration (ng/mL)	O.D 450 nm
0	0.056
5	0.206
15	0.498
50	1.367
100	2.165
200	3.005

This curve was obtained on the date of performance testing of the components. It is for the only purpose of illustration. Each laboratory should establish its own standard curve

CONTROL VALUES

Control No.	Mean (ng/mL)	Range
1	5.89	3.35-8.12
2	21.0	15.8-27.4
3	38.3	30.3-49.0

This lot was tested by test methods described in the manufacturer's package insert.

These results are applicable to this lot and sequential sublots, for which a capital letter have been assigned after the lot number.

This Certificate of Analysis is a computer printout and has therefore not been signed by hand.

This lot meets our specifications and is released for sale	APPROVED	Signed Date of Manufacture, Alícia Abellan Quality Control Dept. 30/08/2017
	(APPROVED / NOT APPROVED)	

