

CERTIFICATE OF ANALYSIS	hCG EIA 96 TEST
--------------------------------	------------------------

REF

6107310	96 Test	X
---------	---------	---

Cont	hCG Ab Coated Wells	1
	Standard Liq	1 ml
	Conjugate reagent	18 ml
	TMB reagent	11 ml
	Stop Solution	11 ml
	Zero Buffer	13 ml

LOT RN-61143 2019-01
 2-8°C **IVD**

STANDARD CURVE

Concentration (mIU/mL)	O.D 450 nm
0	0.045
5	0.126
20	0.385
50	0.888
150	2.093
300	2.999

This curve was obtained on the date of performance testing of the components. It is for the only purpose of illustration. Each laboratory should establish its own standard curve

CONTROL VALUES

Control No.	Mean (mIU/mL)	Range
1	6.06	5.55-10.4
2	17.6	11.1-28.4
3	134	108-167

This lot was tested by test methods described in the manufacturer's package insert.

These results are applicable to this lot and sequential sublots, for which a capital letter have been assigned after the lot number.

This Certificate of Analysis is a computer printout and has therefore not been signed by hand.

This lot meets our specifications and is released for sale	APPROVED	Signed Date of Manufacture, Mercè Oriol Quality Control Dept. 11/09/2017
	(APPROVED / NOT APPROVED)	

