

CERTIFICATE OF ANALYSIS	LH EIA 96 Test
--------------------------------	-----------------------

REF

6107320	96 Test	x
---------	---------	---

Cont	LH Ab Coated Wells	1
	Standard Lyo.	1
	Conjugate Reagent	13 ml
	TMB	11 ml
	Stop Solution	11 ml

LOT RN-56998 2018-10
 2-8°C **IVD**

STANDARD CURVE

Concentration (mIU/mL)	O.D 450 nm
0	0.053
5	0.140
15	0.321
50	0.953
100	1.742
200	2.832

This curve was obtained on the date of performance testing of the components. It is for the only purpose of illustration. Each laboratory should establish its own standard curve

CONTROL VALUES

Control No.	Mean (mIU/mL)	Range
1	2.63	1.13-3.79
2	20.9	16.8-24.5
3	53.5	41.9-63.9

This lot was tested by test methods described in the manufacturer's package insert.

These results are applicable to this lot and sequential sublots, for which a capital letter have been assigned after the lot number.

This Certificate of Analysis is a computer printout and has therefore not been signed by hand.

This lot meets our specifications and is released for sale	APPROVED	Signed Date of Manufacture, Alicia Abellan Quality Control Dept. June 15, 2017
	(APPROVED / NOT APPROVED)	

