

CERTIFICATE OF ANALYSIS	FSH EIA 96 Test
--------------------------------	------------------------

REF	6107325	96 Test		
-----	---------	---------	--	--

Cont	FSH Ab Coated Wells	1
	Standard Lyo.	1
	Conjugate reagent	13 ml
	TMB	11 ml
	Stop Solution	11 ml

LOT	RN-56940	2018-09
	2-8°C	IVD

STANDARD CURVE

Concentration (mIU/mL)	O.D 450 nm
0	0.063
5	0.118
15	0.249
50	0.880
100	1.913
200	3.072

This curve was obtained on the date of performance testing of the components. It is for the only purpose of illustration. Each laboratory should establish its own standard curve

CONTROL VALUES

Control No.	Mean (mIU/mL)	Range
1	9.64	7.22-13.0
2	20.2	16.6-31.6
3	45.4	37.0-73.9

This lot was tested by test methods described in the manufacturer's package insert.

These results are applicable to this lot and sequential sublots, for which a capital letter have been assigned after the lot number.

This Certificate of Analysis is a computer printout and has therefore not been signed by hand.

This lot meets our specifications and is released for sale	APPROVED	Signed Date of Manufacture, Alicia Abellan Quality Control Dept. 15 June 2017
	(APPROVED / NOT APPROVED)	

