

CERTIFICATE OF ANALYSIS	TESTOSTERONE
--------------------------------	---------------------

REF

6107615	96 Test		
---------	---------	--	--

Cont	Goat Anti-Rabbit Ab Coated Wells	1
	Standard set	1
	Controls 1/2	1
	Conjugate reagent	12 ml
	Antibody Reagent	7 ml
	TMB reagent	11 ml
	Stop Solution	11 ml

LOT RN-56747 2018-08
 2-8°C **IVD**

STANDARD CURVE

Concentration (pg/mL)	O.D 450 nm
0	2.506
0.1	2.072
0.5	1.393
2.0	1.002
6.0	0.803
18.0	0.485

This curve was obtained on the date of performance testing of the components. It is for the only purpose of illustration. Each laboratory should establish its own standard curve

CONTROL VALUES

Control No.	Mean (pg/mL)	Range
1	0.563	0.326-0.725
2	10.9	7.89-12.6

This lot was tested by test methods described in the manufacturer's package insert.

These results are applicable to this lot and sequential sublots, for which a capital letter have been assigned after the lot number.

This Certificate of Analysis is a computer printout and has therefore not been signed by hand.

This lot meets our specifications and is released for sale	APPROVED	Signed Date of Manufacture, Alicia Abellan Quality Control Dept. June 15, 2017
	(APPROVED / NOT APPROVED)	

