



<b>CERTIFICATE OF ANALYSIS</b>	<b>Linear URS-10</b>
--------------------------------	----------------------

<b>REF</b>	7101010	100 Test strips
------------	---------	-----------------

<b>LOT</b>	URS7050239	 2019-06
	2-30°C	<b>IVD</b>

### FUNCTIONAL TEST WITH TEST URINE

ANALYSIS		SPECIFICATIONS (Color card)	TEST RESULTS
S.G.	1,000	$\geq 1,000 < 1,005$	Complies
	1,025	$> 1,020$	Complies
pH	5,0	$\geq 5,0 < 6,0$	Complies
	9,0	$> 8,0 \leq 9,0$	Complies
Leukocyte	Neg. Urine	Negative	Complies
	Yihua Urine Controls	$\geq 2+$	Complies
Nitrite	Neg. Urine	Negative	Complies
	Yihua Urine Controls	Positive	Complies
Protein	Neg. Urine	Negative	Complies
	Yihua Urine Controls	$> 1+$	Complies
Glucose	Neg. Urine	Negative	Complies
	Yihua Urine Controls	$\geq 2+$	Complies
Ketone	Neg. Urine	Negative	Complies
	50 mg/dl	$> 1+$	Complies
Urobilinogen	Neg. Urine	Negative	Complies
	4 ~ 12 mg/dl	$> 4$ mg/dl	Complies
Bilirubin	Neg. Urine	Negative	Complies
	Yihua Urine Controls	$> 1+$	Complies
Blood	Neg. Urine	Negative	Complies
	Yihua Urine Controls	$> 2+$ (Hb)	Complies

This lot was tested by test methods described in the manufacturer's package insert.

These results are applicable to this lot and sequential sublots, for which a capital letter have been assigned after the lot number.

This Certificate of Analysis is a computer printout and has therefore not been signed by hand.

<b>This lot meets our specifications and is released for sale</b>	APPROVED	<b>Signed -</b> <b>Date of Manufacture,</b> Alicia Abellan Quality Control Department, 25 July 2017
	(APPROVED / NOT APPROVED)	

